

Suction Cup Grab Bars

HCPCS: E0241



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Installs and removes without tools or professional installation
- Can be installed at any angle
- Large suction cups provide an extremely strong hold suction cups are 3.75" in diameter
- Release levers make installing and removing the suction cup grab bars quick and easy
- Red/Green color indicators on top of the handle shows "safe secure" suction or to "re-affix" suction
- Avoid grout lines for proper placement
- Will not damage property
- Use to assist with balance only
- Do not use to leverage body weight
- Always check to make sure the bar is securely attached before using
- Uninstall and reinstall as needed
- Retail Packaged

PRODUCTS

Description

RTL13083	Grab Bar,Suction,WHT/BEIG,12",1EA,3/CS
RTL13085	Grab Bar,Suction,BLK/CHRM,12",1EA,3/CS

SPECIFICATIONS

Description

Carton Shipping Weight	3.5 lbs.
Length	12"
Warranty	Limited Lifetime
Weight	.5 lbs.
Grab Bars	Length 12"