

## Swivel Seat Shower Stool



INSET A



INSET B



RETAIL BOX

- The rotating padded seat reduces twisting and reaching while in the shower
- Seat rotates 360° and locks every 90° (Figure A)
- The rotating top of the seat can be removed for easy cleaning
- Easy tool-free assembly
- Ideal for individuals who have balance or mobility issues
- Seat height can be adjusted from 16" to 22", in 1" increments
- Retail packaged

### PRODUCTS

### Description

RTL12061M	4 retail boxes/cs
-----------	-------------------

### SPECIFICATIONS

### Description

Outside Legs	15" (W) x 15" (D)
--------------	-------------------

Seat Dimensions	14" (W) x 14" (D)
-----------------	-------------------

Seat Height	16" – 22"
-------------	-----------

Weight	6.6 lbs.
--------	----------

Carton Shipping Weight	28.65 lbs.
------------------------	------------

Weight Limit	450 lbs.
--------------	----------

Warranty	Limited Lifetime
----------	------------------