



Item #	Seat	UOM
TITAN18P22	18' x 18" Full Back Pan Seat	1/cs
TITAN20P22	20' x 20" Full Back Pan Seat	1/cs

accessories		
AA4100	Swing-Away Joystick	1/ea
AE2500	ELR Bracket	1/ea
AH1000	Cane/Crutch Holder	1/ea
AB1000	Arm Bag	1/ea
AH2000	Oxygen Tank Holder	1/ea
AB1120	Backpack	1/ea
AZ0060	Cup Holder	1/ea



HCPCS  
**K0822**



## Titan P22

### Standard Powerchair

- Height-adjustable swivel seat features fold-down backrest.
- Includes interchangeable color panels in red and blue.
- Seat back reclines and features adjustable headrest.
- Adjustable-length controller mount. (Figure B)
- Armrests are padded, adjustable and removable.
- Available with 18" and 20" wide Pan Seats. (Figure A)
- Optional swing-away controller arm (Item #AA4100) sold separately.
- Solid, flip-up footplate with adjustable height, depth and angle.
- Positioning belt included.
- Optional elevating leg rests (requires Item #AE2500 & Item #LK3JELR) sold separately.
- Programmable controller.
- Flat-free, non-marking tires. (Figure C)
- Easy freewheel operation. (Figure D)
- Anti-tip wheels.

#### specifications

HCPCS	K0822
Weight Capacity	300 lbs.
Top Speed	4 mph
Maximum Range	15 miles
Turning Radius	25"
Climbing Angle	6°
Ground Clearance	2.5"
Dimensions	40" (L) x 24" (W)
Floor -to-Seat Height	20.5" - 23"
Controller	Dynamic Shark 40A
Motors (2)	24V x 200W x 3900 RPM
Batteries (pair)	12V x 35AH (U1)

Charger	3A Offboard
Brakes	Electromagnetic
Caster Wheels	6" x 2" Solid
Drive Wheels	10" x 3" Flat-free
Base Weight	68 lbs.
Battery Weight (pair)	44 lbs.
Seat Weight	43 lbs.
Total Weight	158 lbs.
Warranty on Frame	Limited Lifetime
Warranty on Electronic Controller/Drive Train Components	14 months
Warranty on Batteries	6 months