




Item #	Seat	UOM
TITANLB18CSX23	18" x 18" Captain Seat, Standard Back	1/cs
TITAN18CSX23	18" x 18" Captain Seat, Full Back	1/cs
TITAN20CSX23	20" x 18" Captain Seat, Full Back	1/cs
TITAN1616X23	16" x 16" Captain Seat, Full Back	1/cs
TITAN1618X23	16" x 18" Captain Seat, Full Back	1/cs
TITAN1816X23	18" x 16" Captain Seat, Full Back	1/cs
TITAN2020X23	20" x 20" Captain Seat, Full Back	1/cs
TITAN22CSX23	22" x 20" Captain Seat, Full Back	1/cs



HCPCS
K0823

accessories		
AA4100	Swing-Away Joystick	1/ea
AE2500	ELR Bracket	1/ea
LK3JELR	Plastic Foot plate	1/ea
AH1000	Cane/Crutch Holder	1/ea
AH2000	Oxygen Tank Holder	1/ea
AB1000	Arm Bag	1/ea
AB1120	Backpack	1/ea
AZ0060	Cup Holder	1/ea
ST205-COVER	18" Seat Cover	1/ea
ST301-COVER	20" Seat Cover	1/ea
ST306-COVER	22" Seat Cover	1/ea



Titan X23

Standard Powerchair

- Height-adjustable swivel seat features fold-down backrest. (Figure A)
- Includes interchangeable color panels in red and blue.
- Seat back reclines and features adjustable headrest.
- Adjustable-length controller mount. (Figure B)
- Armrests are padded, adjustable and removable.
- Available with 16", 18", 20" and 22" wide Captain's Seat.
- Optional swing-away controller arm (Item #AA4100) sold separately.
- Solid, flip-up footplate with adjustable height, depth and angle
- Positioning belt included.
- Optional elevating leg rests (requires Item #AE2500 & Item #LK3JELR) sold separately.
- Programmable controller.
- Flat-free, non-marking tires. (Figure C)
- Easy freewheel operation. (Figure D)
- Anti-tip wheels

specifications

HCPCS	K0823
Weight Capacity	300 lbs.
Top Speed	4 mph
Maximum Range	15 miles
Turning Radius	25"
Climbing Angle	6°
Ground Clearance	2.5"
Dimensions	40" (L) x 24" (W)
Floor -to-Seat Height	20.5" - 23"
Controller	Dynamic Shark 40A
Motors	(2) 24V x 200W x 3900 RPM
Batteries (pair)	12V x 35AH (U1)

Charger	3A Offboard
Brakes	Electromagnetic
Caster Wheels	6" x 2" Solid
Drive Wheels	10" x 3" Flat-free
Base Weight	68 lbs.
Battery Weight (pair)	44 lbs.
Seat Weight	43 lbs.
Heaviest Piece Disassembled	55 lbs.
Total Weight	158 lbs.
Warranty on Frame	Limited Lifetime
Warranty on Electronic Controller/Drive Train Components	14 months
Warranty on Batteries	6 months