



## Toilet Safety Frame with Padded Arms



FIGURE A



FIGURE B



RETAIL BOX

- Anodized aluminum is sturdy and lightweight
- Padded armrests provide additional comfort and support
- Comes with height-adjustable (1" increments) and width-adjustable arms (2 positions) (Figure A)
- Powder-coated aluminum bracket easily attaches frame to toilet (Figure B)
- Bracket secures safety frame in place

### PRODUCTS

PRODUCTS	Description
RTL12000	Unassembled, 1/cs

### SPECIFICATIONS

SPECIFICATIONS	Description
Width Between Arms	19" (smallest)   23.50" (largest)
Adjustable Height	28" (smallest)   32" (highest)
Weight	4 lbs.
Carton Shipping Weight	20 lbs.
Weight Limit	300 lbs.
Depth of Frame	18.5"
Warranty	Limited Lifetime