

## Toilet Safety Rail



RETAIL BOX

- Allows individuals to safely sit or stand
- Easy, tool-free assembly and installation
- Comes with soft, padded handles for comfort
- Width adjusts from 16.5" - 19.5" between handles

### PRODUCTS

### Description

RTL12087	2 retail boxes/cs
----------	-------------------

### SPECIFICATIONS

### Description

Width Between Arms	16.5" - 19.5"
Weight	10.5 lbs.
Carton Shipping Weight	24 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime