

Premium Raised Toilet Seat with Removable Arms



RETAIL BOX

- Designed for individuals who have difficulty sitting down or getting up from the toilet
- Allows individual to use existing toilet seat and lid
- Tool-free arms are 19.75" apart and can be removed or added as needed
- Arm sleeves and receivers are made of metal
- Heavy-duty molded plastic construction
- Lightweight and portable
- Easily attaches to toilet bowl
- No tools required for installation
- Easy to clean
- Two models accommodate standard and elongated toilets

PRODUCTS

Description

12402	Fits Standard Toilets, 1/cs
12403	Fits Elongated Toilets, 1/cs

SPECIFICATIONS

Description

Seat Dimensions (12402)	13.75" (W) x 17" (D) x 3.5" (H)
Seat Dimensions (12403)	13.75" (W) x 19.25" (D) x 3.5" (H)
Width Between Arms	19.75"
Weight	3.5 lbs.
Carton Shipping Weight	4.5 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime