



Item #	Seat	UOM
2850-18	18" Captain's Seat	1/cs
2850-20	20" Captain's Seat	1/cs



HCPCS
K0823

Trident

Standard Power Wheelchair Front-Wheel Drive

- Shroudless design offers easy access to internal components and batteries.
- Semi-reclining seat with adjustable headrest offers three height adjustments and two depth options.
- Armrests are padded, adjustable and removable.
- Available with 18" or 20" width Captain's Seat.
- Quick-release pins allow easy seat removal.
- Adjustable, flip-up footplate includes rubber front bumper to prevent damage.
- Positioning belt included.
- Seat position can be adjusted forward and backward.
- Programmable controller.
- Optional swing-away controller arm sold separately.
- Optional elevating leg rests sold separately.
- Black, flat-free non-marking wheels.
- Large, front anti-tip wheels feature suspension for safety, smoother ride and enhanced obstacle handling.
- Controller is protected by metal frame.
- Easy freewheel operation.
- Slide-out battery tray.

specifications

HCPCS	K0823	Charger	4A Offboard
Weight Capacity	300 lbs.	Brakes	Electromagnetic
Top Speed	4 mph	Caster Wheels	8" x 2" Solid
Maximum Range	15 miles	Drive Wheels	14" x 3" Flat-free
Turning Radius	24"	Base Weight	82 lbs.
Climbing Angle	7.5°	Battery Weight (pair)	53 lbs.
Ground Clearance	2"	Seat Weight	51 lbs.
Dimensions	42.25" (L) x 27.25" (W)	Total Weight	186 lbs.
Floor-to-Seat Height	21" - 24"	Seat Dimensions	18" (W) x 18" (D) 20" (W) x 18" (D)
Controller	PG VR2 50A	Warranty of Frame	Lifetime
Motors (2)	24V x 320W x 4600 rpm	Warranty on Electronic Controller/Drive Train Components	14 Months
Batteries (pair)	12V x 36AH (U1)	Warranty on Batteries	6 months