



### Deluxe, Trigger Release Folding Walker with 5" Wheels



FIGURE A

- Trigger release feature for opening and closing allows users to close the walker without removing their hands from the handle and is ideal for individuals with limited finger dexterity
- Trigger release can be reached with hands of all sizes Patented
- Composite Trigger Release resists cracking
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength
- Newly designed rear glide cap allows for walker to glide easily and smoothly over most surfaces (Figure A)

#### PRODUCTS

#### Description

10226-4	Adult, Flame Blue Crossbrace, 4/cs
10227-4	Junior, Flame Blue Crossbrace, 4/cs

#### SPECIFICATIONS

#### Description

Width	18" inside hand grip   20.75" inside back legs
Overall Width	25"
Depth (Opened at Base) (10226   10227)	16.75"   16.25"
Weight (10226   10227)	7/6 lbs.
Carton Shipping Weight (10226   10227)	29   27 lbs.
Weight Limit	350 lbs.
Warranty	Limited Lifetime