

## Universal 5" Walker Wheels Dual-Hole Adjustment Column and Two Sets of Rear Glides



RETAIL BOX

- Converts folding walker into wheeled walker
- Fits most manufacturers' walkers
- Allows for 8 height adjustments
- Rubber wheel allows walker to roll easily and smoothly over irregular surfaces
- Comes with rear glide caps, (item #10107) and glide covers (item #10107C) allowing use on all surfaces
- 350 lb Weight Limit

### PRODUCTS

	Description
10109	Silver tubing, gray tire and silencer, 1/pr

### SPECIFICATIONS

	Description
Weight Limit	350 lbs.