



COMPETITIVE
EDGE Line



Durable 4 Wheel Rollator with 7.5" Casters

Fold-Up and Removable Back Support, Loop Locks

Item #	Color	UOM
10257BL-2	Blue	2/cs
10257RD-1	Red	1/cs

- Comes with durable, padded seat.
- Removable, hinged, padded backrest can be folded up and down.
- Comes with a basket to store personal items.
- Durable steel frame.
- 7.5" non-marring casters are ideal for indoor and outdoor use.
- Easy-to-use deluxe loop locks. (Figure A)
- Lever allows for independent handle height and angle adjustment. (Figure B)

specifications

Handle Height	31" - 37"
Unit Dimensions	23.5" (W) x 25.5" (L)
Seat Dimensions	15" (W) x 21" (H) x 8.5" (D)
Weight	20 lbs.
Carton Shipping Weight	26 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

accessories

809 Tray	1/ea
----------	------

Aluminum Rollator, 6" Casters

Fold-Up and Removable Back Support, Padded Seat, Loop Locks

Item #	Color	UOM	Item #	Color	UOM
R726BL	Blue	1/cs	R726GR	Green	1/cs
R726RD	Red	1/cs	R726BK	Black	1/cs
R726SL	Silver	1/cs			

- Comes with new seamless padded seat (Figure A).
- Zippered pouch under seat for added privacy and security of personal items. (Figure B)
- Removable, hinged, padded backrest can be folded up or down.
- 6" Black non-marring casters with soft grip tires are ideal for indoor and outdoor use. Easy-to-use deluxe loop locks. (Figure C)
- Brakes with serrated edges provide firm hold. (Figure D)

specifications

Handle Height	32" - 37"
Unit Dimensions	24.5" (W) x 27.5" (L)
Seat Dimensions	12" (W) x 20" (H) x 14" (D)
Weight	13 lbs.
Carton Shipping Weight	16 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

accessories

9501S72606 Basket	1/ea
-------------------	------