

ZooMe 3-Wheel Recreational Scooter



HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

The ZooMe 3 is a modern recreational electric scooter providing personal transportation for everyone! The comfortable extra-wide seat is removable, so you can stand up or sit down while riding it making it great for recreational or mobility use. The ZooMe 3 comes with a powerful 350 watt brushless hub motor providing a top speed of 15 mph, and a maximum range of 17 miles on a full charge.

- Powerful hub driven motor providing speeds up to 15 mph
- Folding tiller for easy transportation in a vehicle, or for storage
- Two keys included
- Hand brake with locking mechanism for parking on inclines
- Forward and reverse controls are easily changed with the switch of a button
- Fast and slow speed settings easily changed with the switch of a button
- Front basket for storage
- Pneumatic tires for the smoothest ride possible
- Extra bright LED headlight
- Anti-tip wheel

PRODUCTS

PRODUCTS	Description
ZOOME3	ZooMe 3-Wheel Rec. Scooter, Red, 1/ea

Description

SPECIFICATIONS

SPECIFICATIONS	Description
Warranty	24 Month Limited Warranty on Electronics
Product Weight	78 lb
Product Dimensions	41" (L) x 27" (W) x 29" (H)
Seat Dimensions	10" (D) x 14.5" (W)
Weight Capacity	300 lb
Brakes	Hand Brake
Top Speed	15 mph
Anti-Tip Wheels	Standard
Batteries (pair)	12V x 12AH (3)
Battery Charger	Offboard 1.6A Smart Charger
Climbing Angle	10°
Color	Red
Controller	N/A
Ground Clearance	4.25"
HCPCS	N/A
Maximum Range	20 miles
Motor	350W Brushless Hub

SPECIFICATIONS**Description**

Seat Type	Removable - Post Style
Turning Radius	22"
Warranty on Batteries	12 Month Warranty on Batteries
Warranty on Frame	Limited Lifetime Warranty on Frame
Wheel Size	Front: 16" x 2" Rear 4" x 2"