

ZooMe Auto-Flex Folding Travel Scooter



HCPCS:



FIGURE A



FIGURE B




FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

The Drive Medical ZooMe Auto-Flex is ideal for on-the-go lifestyles combining a lightweight design and premium "Flex" feature of single-motion folding and unfolding for easy transportation. Unlike most travel scooters, no assembly or disassembly required! The compact frame allows you to easily transport or store the scooter in the trunk of a car. The 4-wheel design provides optimal stability on all types of terrain. The ZooMe Auto-Flex electric scooter includes an angle-adjustable delta tiller for maximum driving comfort, and ensures easy throttle control for users with limited dexterity.

- Fully electronic folding and unfolding in 15 seconds or less with the wireless key fob
- Aluminum alloy frame is lightweight for superior portability
- Airline safe lithium battery
- High-gloss shrouds for a protective finish
- Easy throttle control
- Delta tiller for users with limited dexterity
- Auto-folding seat with padded, adjustable armrests
- Flat-free, non-marking tires for worry-free travel
- Easy-to-adjust tiller places the scooter's controls in the ideal position
- Secondary manual-release levers and toggle switch for folding and unfolding
- Comes with two ignition keys and one key fob included
- Front and rear reflectors
- Anti-tip wheels for additional safety

PRODUCTS	Description
FLEX-AUTO	ZooMe Flex Auto Folding Scooter, Blue
SPECIFICATIONS	Description
Warranty	24 Month Limited Warranty on Electronics
Product Weight	60 lb
Product Dimensions	37" (L) x 17" (W) x 33" (H)
Seat Dimensions	14" (D) x 16" (W)
Weight Capacity	300 lb
Brakes	Electromagnetic
Top Speed	4 mph
Anti-Tip Wheels	2" x .5"
Batteries (pair)	25.9V x 10AH Lithium
Battery Charger	2A Offboard
Battery Weight (pair)	5 lb
Climbing Angle	6°
Color	Blue

SPECIFICATIONS	Description
Controller	PG S Drive 45A
Freewheel Mode	Yes
Ground Clearance	1.6"
HCPCS	N/A
Heaviest Piece Weight	55 lb
Maximum Range	13 miles
Motor	24V x 120W
Seat Type	Stadium
Turning Radius	47.2"
Warranty on Batteries	12 Month Warranty on Batteries
Warranty on Frame	Limited Lifetime Warranty on Frame
Wheel Size	Front: 6" x 1.5" Rear: 7" x 2.35"