

## Silver Sport Recliner with Vinyl Upholstery

HCPCS:



- Vinyl upholstery is comfortable, lightweight and easy-to-clean
- Hydraulic reclining mechanism allows for infinite adjustments up to 180°
- Arms include support extension for reclining position
- Carbon Steel frame with silver vein finish is attractive, and easy-to-maintain
- Composite Mag-style wheels are lightweight and maintenance free
- 8" front casters
- Padded armrests provide added patient comfort
- Wheels are set back on frame to prevent tipping
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Rear anti-tippers come standard
- Headrest extension with cushioned head immobilizer come standard
- Includes swing-away elevating leg rests with Metal foot plates
- Comes with push to lock wheel locks

PRODUCTS	Description
SSP16RBDDAV	SilverSport16"ReclBack,DetDesk,Vnl,1/ea
SSP16RBDFAV	SilverSport16"ReclBack,FullArm,Vnl,1/ea
SSP18RBDDAV	SilverSport18"ReclBack,DetDesk,Vnl,1/ea
SSP18RBDFAV	SilverSport18"ReclBack,FullArm,Vnl,1/ea
SSP20RBDDAV	SilverSport20"ReclBack,DetDesk,Vnl,1/ea
SSP20RBDFAV	SilverSport20"ReclBack,FullArm,Vnl,1/ea

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	59 lb   60 lb   61 lb
Product Dimensions	50" (L) w/ riggings x 23" (W) x 50" (H)
Material	Steel frame   Vinyl upholstery
Seat Dimensions	16" (W) x 16" (D)
Weight Capacity	16" - 18" 300 lb   20" 350 lb
Seat To Floor Height	20.5"
Carton Shipping Weight	66 lb   67 lb   67.5 lb
Folded Dimensions	50" (L) x 12" (W) x 50" (H)
CA Prop 65	Yes
Back Height	21"   31" w/ headrest

ACCESSORIES	Description
STDS856	Seat Belt,60"Velcro-type,1/ea, EA

ACCESSORIES	Description
STDS850	Seat Belt, Clasp Type,1/ea, EA
STDS855	Seat Belt,60" Auto-type,1/ea, EA
STDS1040S	Cup Holder, Swivel,1BX,6/CS, BX
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
STDS818	Anti-Tippers,rear,SentraRecliner,1pr/bx, BX
STDELR-AL	ArticulatingElev.Legrest, Std,1pr/cs, CS
LK3JELR	ELR Set,Black,1Pair eq EA, PR
WASR	Limb Support, Right,1/cs, CS