

NEW!

COMPARE TO:
COMPARE TO INVACARE® SOLACE & SPAN AMERICA GEOMAT HC



15019


Available 1 Pack Therapeutic 5 Zone Support Mattress

15019 80" (L) x 36" (W) x 6" (H), 1/ea
HCPCS code: E0184

FEATURES:

- Deluxe, horizontal, cross-cut, **5-zoned foam** mattress provides comfort, support and pressure redistribution.
- 210 denier nylon cover reduces friction and shear, and is water resistant and vapor permeable.
- 275 lb. weight capacity.
- Limited 3 Year Warranty on foam.
- Limited 1 Year Warranty on cover.

NEW!



15770

Available 1 Pack Gravity 7 Long Term Care Pressure Redistribution Mattress

15770 80" (L) x 36" (W) x 6" (H), 1/ea
 15777 80" (L) x 36" (W) x 6" (H) with 3" Elevated Perimeter and Cut Out, 1/ea
 15784 84" (L) x 36" (W) x 6" (H), 1/ea
 15785 84" (L) x 36" (W) x 6" (H) with 3" Elevated Perimeter and Cut Out, 1/ea
 15876 76" (L) x 36" (W) x 6" (H), 1/ea
 15877 76" (L) x 36" (W) x 6" (H), with 3" Elevated Perimeter and Cut Out, 1/ea

HCPCS code: E0184

FEATURES:

- Multi-layered, high density foam with die-cut top layer provides optimal pressure redistribution.
- The **ultra-soft heel section** is uniquely designed to provide extremely low skin interface pressure.
- The mattress cover is a **2-way stretch**, water resistant, vapor permeable low shear material.
- **Optional firm side rails** help the patient transfer in and out of bed minimizing the risk of falling.
- Meets fire code 1633.
- 300 lb. weight capacity.
- **Limited 7 Year Warranty.**

SPECIFICATIONS

Item #	Length	Width	Height	Product Weight	Carton Shipping Weight	Weight Limit
15019	80"	36"	6"	30 lbs.	35 lbs.	275 lbs.
15770	80"	36"	6"	30 lbs.	35 lbs.	300 lbs.
15777	80"	36"	6" + 3" EP	30 lbs.	35 lbs.	300 lbs.
15784	84"	36"	6"	18.3 lbs.	23 lbs.	300 lbs.
15785	84"	36"	6" + 3" EP	19.3 lbs.	24 lbs.	300 lbs.
15876	76"	36"	6"	19.3 lbs.	24 lbs.	300 lbs.
15877	76"	36"	6" + 3" EP	17.5 lbs.	22 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.