


This stool is sometimes used as a hip chair, but is more commonly found in bars and lounges...



Generic Hip Chair

ASCENDER HIP CHAIR

...this chair is only found on orthopedic wards



ASCENDER 300
Ultra Hip with
Leg Support

ORTHOPEDIC (HIP) CHAIR COMPARISON

	Generic Hip Chair	ASCENDER Hip Chair	Advantage
Seat Height	27 inches	23 or 25 inches	ASCENDER
Seat Width	19 inches	20, 24, 28 inches	ASCENDER
Seat Depth	17 inches	19 inches	ASCENDER
Leg Support Option	No	Yes	ASCENDER
Slant Seat Option	No	Yes	ASCENDER
Retractable Foot Rest	No	Yes	ASCENDER
A-Frame Design	No	Yes	ASCENDER
Adjustable Back	No	Yes	ASCENDER
Weight Capacity	usually < 300 lbs	up to 1,000 lbs	ASCENDER



ASCENDER
M.P. 500



ASCENDER 500
Transporter

AscenderTM
ORTHOPEDIC AND BARIATRIC
Seating Solutions

Ascender™

ORTHOPEDIC AND BARIATRIC
Seating Solutions

PRODUCT SPECIFICATIONS

Description	Features	Adj. Seat Hgt (in.)	300 lbs 20" W	500 lbs 24" W	1000 lbs 28" W
			Model #	Model #	Model #
Hip Chair	Standard	25/23	7000	7300	7500
Ultra Hip	Leg up Support	25/23	7001	7301	7501
Slant Hip	Forward Sloping Seat	25/23	7100	7700	7800
Ultra Slant Hip	Forward Sloping Seat and Leg up Support	25/23	7101	7701	7801
Transporter	Wheels	25/23	7200	7400	7600
Ultra Transporter	Wheels and Leg Support	25/23	7201	7401	7601
Multi-Purpose Chair	Forward Sloping Seat	20/18	8000	8200	8400
Home Chair	Adjustable Slant	20	9000	9200 (24" wide)	



Multi-Purpose Chair



Hip Chair



Transporter



Ultra Hip



Home Chair

CUSHION COLORS



Wine Rose



Ice Blue



Emerald



Wheat

STAIN COLORS



Light



Dark

ASCENDER and GLIDE-A-WAY are trademarks of Rehab Seating Systems, Inc.
US Patent Nos. 4358156 and 5277477
Patent Pending on M.P. models
Made in USA



Eight Alton Place • Suite 3
Brookline, MA 02446
Phone: 617-734-8500 or 800-525-SEAT (7328)
Fax 617-738-9894
e-mail: info@rehabseating.com
www.rehabseating.com