



# Molift Smart 150

Light, foldable and transportable without tools



## Compact smartness

Molift Smart 150 is designed to be easily wheeled and transported to where it is required. It is designed for home care, for travelling and easy storage. The Molift Smart is foldable in one piece without using any tools.

Due to its low weight, the small footprint and the unique construction, Molift Smart 150 is easy to manoeuvre even in narrow spaces. The hoisting range of 270–1680 mm provides an excellent maximum hoisting height and at the same time gets very low, facilitating hoisting from the floor. The 4-point sling bar ensures a comfortable and spacious hoisting position for the user.

Molift Smart 150 has been awarded for Design Excellence by the Norwegian Design Council.



Product name	Item no.	Dimensions		
Molift Smart 150	1020	Length	1200 mm	
		Leg height	110 mm	
		Dimension folded	L1160 x W475 x H360 mm	
		Hoisting height	270–1665 mm	
		Hoisting interval	1410 mm	
<b>Material</b>		<b>Technical</b>		 <p>Clean / disinfect the product using a standard non-abrasive cleaning agent with a pH between 5 and 9, or with a 70% disinfectant solution. The product can be safely cleaned in a cabinet washer-disinfector at 85°C for 3 minutes.</p>
Aluminium, steel		Service Software	Included	
<b>Weight</b>		Battery	NiMH 14.4 V - 2.6 Ah	
Total Weight	26 kg	Battery charge time	approx. 3 hours	
Weight when separated	Part 1: 13 kg Part 2: 13 kg	Hoisting speed	60 mm/sec	
		Rated performance	50 hoists/75 kg/500 mm	

## Accessories



### Molift Soft Travel Bag

The Molift Soft Travel Bag protects and enables easy transporting of the hoist, battery, charger and sling during travel.

Item no: 3049177



### Travel Suitcase

This hard shell suitcase provides additional protection when transporting by air. It is TSA approved for air travel and ideal for the client with an active life.

Item no: 0990200

## Features



### 4-point sling bar

The standard 4-point sling bar ensures a comfortable and spacious hoisting position for the user.



### Foldable without tools

Molift Smart 150 is designed to be easy and convenient at home and to transport when travelling. The possibility to easily fold the hoist without tools is a big advantage.



For the latest news and continuously updated product information – visit: [www.etac.com](http://www.etac.com)