

All of the features and benefits of the Preferred Care Recliner plus a unique drop arm mechanism. The self-storing, height adjustable drop arm is simple for the caregiver to activate and enables a safer, easier lateral transfer for the resident and caregiver.

Available Colors:

CA TB 133 Upholstery

- Blue Ridge: FR565DG427
- Jade: FR565DG857
- Rosewood: FR565DG863



Lumex Preferred Care Recliner with Tilt-in-Space Function: Series FR565TG

The Tilt-in-Space feature allows the caregiver to change the “tilt” of the seat and back for pressure relief and position change. Please note that this version of the recliner does not recline but utilizes a side mounted foot pedal mechanism to extend the legrest and “tilt” the position of the patient while maintaining a constant back to seat angle. In addition, the recliner features a self-storing, height adjustable drop arm that is simple for the caregiver to activate and enables a safer, easier lateral transfer for the caregiver and resident.

Available Colors:

CA TB 133 Upholstery

- Blue Ridge: FR565TG427
- Rosewood: FR565TG863

