

Standard features found on all Lumex Preferred Care Recliners:

- Reduced armrest and seat height enable the recliner to fit under dining / activity tables and allow residents to “scoot” more easily.
- Two swivel casters with two rear locking casters enable the Preferred Care Recliner to be easily maneuvered in small spaces and track straight when pushed down hallways.
- A unique Ultra-Flex suspension system in the seat and back adjusts to each person, cradling the resident from head to heel, greatly reducing the risk of “bottoming out”.
- By simply removing the back upholstery from the recliner back, the Ultra-Flex Suspension Bands can be adjusted to meet individual patient comfort requirements.
- Extra-wide padded armrests and side panels maximize comfort.
- Quick-release upholstery allows for easy cleaning, removal and replacement.
- A variety of accessories is available to customize the recliners for specific resident needs.
- CA TB-133 Rated Upholstery
- Weight Capacity: 250 lb Evenly Distributed

**Lumex Preferred Care Recliner-Series FR565G****Available Colors:**

CA TB-133 Upholstery

- Blue Ridge: FR565G427
- Jade: FR565G857
- Rosewood: FR565G863