



INVACARE® MICROAIR® MA500  
THERAPEUTIC SUPPORT SURFACE



Yes, you can.®

microAIR®

# Invacare® microAIR® MA500

Alternating Pressure Low Air Loss Mattress System

**Low - Medium Risk  
Stage I-III**

A lightweight and easy-to-use Alternating Pressure Mattress System with Low Air Loss aids in the Prevention and Treatment of Pressure Injuries



**Highly vapor permeable and oversized pliable quilted nylon top cover** provides low shear, friction and moisture protection



**Cable management system** keeps cords away from the ground to help prevent trip hazards




**Nylon non-slip base cover** helps reduce shear and friction during bed articulation



# Features & Benefits

**8 adjustable comfort settings** depending on the patient's therapy requirement



ALTER  
STATIC

- POWER
- POWER FAILURE
- LOW PRESSURE
- STATIC

ON

OFF

**microAIR 500**

**2-in-1 alternating pressure cycles** help ensure protection for patients with low to medium risk



**Soft nylon urethane backed cells** provide easy cleaning and helps minimize the proliferation of bacterial growth

**Cell in cell design** helps prevent bottoming out in the event of a power failure

**Multiple snap bladders** make cells more secure with patient movement

**Natural heel slope** helps decrease pressure in patient's vulnerable heel area

# Request a Quote or Place an Order Today!



**Contact your Inside Sales Representative**

or



**Contact Customer Service at 800.333.6900**

**Invacare Corporation**  
www.invacare.com

**USA**  
One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

**Canada**  
570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

## microAIR MA500 Alternating Pressure Low Air Loss Mattress System



Model #	Dimensions	Weight Capacity	Flow Rate (Pump)	Therapy Type	Alarms	Limited Warranty
MA500	80" L 36" W 8" H	350 lb.	5 LPM	Alternating Pressure Low Air Loss	Audio & Visual	1 Year (Cover/Mattress) 1 Year (Power Unit)
MA500P	microAIR MA500 Pump only					
MA500M	microAIR MA500 Mattress only					

©2019  
Invacare Corporation.  
All rights reserved.  
Trademarks are identified by the symbols™ and ©. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.  
Form No. 19-022  
190201

