



Invacare® Premier Series  
Roze™ Stand-Up Lift

---

***Yes, you can.®***



We leave the lifting to the lift. For residents who need assistance in standing, Roze stand-up model lifts are expressly designed to promote independence for the resident while at the same time relieving the lifting burden from the operator.

# Invacare® Premier Series Roze™ Stand-Up Lift

---



## Strength Capacity

The unique height-adjustable kneepad and 450 pound weight capacity easily accommodates the majority of residents in your facility.



## Accommodating

For those who only need help rising, simply remove the footplate.



## Ease of Use

With a large push handle, lighter steel and aluminum construction, and simple sling attachment, it's easy to accomplish everyday transfers.



## Maneuverability

A small footprint combined with lighter weight aluminum and steel construction provide easy maneuverability.

## Roze™ Specifications

- Height at Sling Hook-up - MAX : 63.5 inches
- Height at Sling Hook-up - MIN : 37 inches
- Base Width OPEN : 44 inches
- Base Width CLOSED : 26 inches
- Base Height (Clearance) : 4.8 inches
- Base Length : 38.8 inches
- Overall Height : 49.1 inches
- Overall Length : 42 inches
- Overall Width : 26 inches
- Caster Size (Front/Rear) : 4.0 inches
- Weight Capacity : 450 lbs/204 Kg
- Weight Out of Carton : 117 lbs
- Battery : 24V DC (Rechargeable) 2 Shipped with Unit
- Charger Input : 100-240V AC
- Charger Output/Charging Time : 29.5 DC Max 6 hours
- Audio/Visual Low Battery Alarm : Yes
- Motor Safety Devices : Anti-Entrapment
- Limited Warranty : 3 Years On Lift, 2 Years on Electronics
- Emergency Stop Button : Yes
- Emergency Up & Down : Yes



**Warning:** Before using the Patient Sling and Patient Lift, read and understand the Owner's Manual for proper operation and safety procedures.

**Invacare Continuing Care, Inc.**  
A Subsidiary of Invacare Corporation

[www.invacare-ccg.com](http://www.invacare-ccg.com)

[info\\_iccg@invacare.com](mailto:info_iccg@invacare.com)

### USA

One Invacare Way  
Elyria, Ohio 44035  
USA  
800-668-2337

### Canada

994 Hargrieve Road  
London, Ontario  
N6E 1P5 Canada  
800-668-2337

Trademarks are identified by the symbol™ and®.  
All trademarks are owned or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.



**Yes, you can.®**