

INVACARE®  
**TRACER® SX5**  
SERIES WHEELCHAIRS



The **Invacare® Tracer® SX5 Wheelchair's** lightweight frame, weighing less than 36 lb., makes this chair perfect for rental and long- or short-term use. This wheelchair offers multiple seat widths, arm styles, sectional back and seat-to-floor heights by utilizing the dual axle position frame base. The 20" and 22" widths include a semi-pneumatic caster to accommodate a weight capacity of 300 lb. Combining the design and technology of the Invacare® 9000 series wheelchairs make it possible to interchange components with the Tracer SX5 wheelchair.

The **Invacare® Tracer® SX5 Recliner Wheelchair** offers the durability of a lightweight manual wheelchair frame with the added versatility and comfort of a full reclining back. The Tracer SX5 Recliner wheelchair offers dynamic recline range from 90° to 180°. This chair is available with multiple seat widths, arm styles, front riggings, a full range of accessories and a 300 lb. weight capacity on 20" and 22" widths.

### Invacare® Tracer® SX5 Wheelchair Features

Durable, low-maintenance, carbon steel frame is long-lasting  
 Dual axle positions allow variation of seat-to-floor heights  
 14-gauge crossbraces add strength and durability  
 Urethane rear tires mounted on "no flex" wheels offer superior performance  
 Sectional back available in 16" and 18" back heights  
 Low-maintenance, dual-sealed precision bearings  
 Heavy-duty inner liners keep seat and back from stretching  
 Meets RESNA standards for manual wheelchairs in product performance

### Invacare® Tracer® SX5 Recliner Wheelchair Features

Dynamic range of recline from 90° to 180°  
 Durable, low-maintenance, carbon steel frame  
 Urethane rear tires, mounted on "no flex" wheels, offer superior performance  
 Low-maintenance, dual-sealed precision bearings  
 Heavy-duty inner liners keep seat and back from stretching  
 Upholstery is dual-embossed  
 Meets RESNA standards for manual wheelchairs in product performance



1. Dual axle positions
2. Urethane rear tires
3. Headrest
4. Seating options available
5. Builder options available



**Invacare Corporation**  
 www.invacare.com

#### USA

One Invacare Way  
 Elyria, Ohio  
 44035-4190  
 (800) 333-8800

#### Canada

570 Matheson Blvd. E.,  
 Unit 8  
 Mississauga, Ontario  
 L4Z 4G4 Canada  
 (800) 668-5324

©2019

Invacare Corporation.  
 All rights reserved.  
 Trademarks are identified by the symbols ™ and ©.  
 All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.  
 Form No. 11-225  
 Rev. 04/14  
 131066

Model #	Overall Width	Overall Height	Overall Depth w/o riggings	Folded Width	Seat Dimensions	Seat-to-Floor Height	Back Height
TRSX50, TRSX52 TRSX56, TRSX58	22.5"-30.5"	N/A	34"	N/A	14"-22"x16"	17.5", 19.5"	16", 18"
TRSX5RC6, TRSX5RC8	N/A	N/A	30.25"	N/A	14"-22"x16"	17.5", 19.5"	34"

  

	Caster Height	Wheel Height	Weight w/o Rigging	Weight Capacity	Shipping Weight	Carton Dimensions	Limited Warranty
TRSX50, TRSX52 TRSX56, TRSX58	6", 8"	20", 22", 24"	34 lb.	250 lb./300 lb.	43 lb.	32"x13"x 37"	5 Years
TRSX5RC6, TRSX5RC8	8"	20", 22", 24"	53 lb.	250 lb./300 lb.	62 lb.	36"x13"x 44"	5 Years



**Yes, you can.™**