

ProBasics® Knee Walker

HCPCS CODE: E0118

The perfect aid for those recovering from ankle or foot conditions such as surgery, breaks, sprains and amputations, the ProBasics Knee Walker is comfortable to use and easy to maneuver. The handle and seat height can be easily adjusted without tools to meet a variety of user needs. The Knee Walker is simple to fold for storage and comes with a removable wire storage basket.

Model #	Overall Dimensions		Kneepad Dimensions	Kneepad-to-Floor Height	Handle Height Adjustment	Caster Dimensions	Weight Capacity	Units Per Carton	Limited Warranty
1030KW	16"W 31"L	6"W 14"L	17"-23"	32"-41.75"	7.5"	300 lb.	1	3 Years Frame 1 Year Seat and Brakes	

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

