

ProBasics® Bariatric Bed

ProBasics® Full-Electric Bariatric Bed

- Heavy steel plate with pan deck design ensures extra strength and user safety
- Hand unit electrically controls positioning for upper body, knees and bed height to provide maximum comfort and flexibility for both the user and caregiver
- HCPCS Code: E0301

Model #	Description
PBHB4	Full-Electric Bariatric Bed



ProBasics® Bed Rails

A. ProBasics™ Four Bar Half Rails

- Clamp on rail, no tools required
- Welded steel construction with brown vein finish and chrome on the second and fifth rail for scratch-resistant movement of the adjustment bracket
- Compatible with VC5210*, VC5310, 5307IVC, 5310IVC and 5410IVC
- HCPCS Code: E0305
- Invacare and ProBasics beds, mattresses and rails, excluding powered support surfaces and bariatric beds, when used as a system and sold new and unused after August 1, 2007, meet the voluntary FDA Guidelines for bed rail entrapment. Visit www.invacare.com/bedsafety for additional safety and entrapment information.

* Only compatible when used with casters

Model #	Description
PB7035	Four Bar Half Rail



B. ProBasics® Reduced Gap Full-Length Rails

- Welded steel construction with brown vein finish
- Self adjusting spring loaded crossbraces
- Four crossbars help reduce entrapment hazards
- Compatible with VC5210, VC5310, 5307IVC, 5310IVC and 5410IVC
- HCPCS Code: E0310
- Invacare and ProBasics beds, mattresses and rails, excluding powered support surfaces and bariatric beds, when used as a system and sold new and unused after August 1, 2007, meet the voluntary FDA Guidelines for bed rail entrapment. Visit www.invacare.com/bedsafety for additional safety and entrapment information.

Model #	Description
PB6034	Reduced Gap Full-Length Rail



Model #	Dimensions	Sleep Surface	Weight Capacity	Shipping Weight	Carton Dimensions	Limited Warranty
PBHB4	88"L	42"W x 80"D	600 lb.	N/A	N/A	2 Years Frame
PB3000 (headspring)	18.25" x 26.25"H			104 lb.	45" x 10" x 43.5"	1 Year Electronics/Motors
PB3010 (footspring)	42"W			138 lb.	47" x 10" x 44.5"	
PB3301 (bed ends)				71 lb.	45" x 9.25" x 23.5"	
PB7035	37"L x 19"H	N/A	N/A	23 lb.	37" x 7" x 17"	2 Years
PB6034	55"L 13.5"H Raised 4.5"H Lowered	N/A	N/A	27 lb.	41.5" x 3" x 19"	2 Years

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.