



Yes, you can.®

Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

Invacare® Solara™ 3G - Adult Tilt-in-Space Wheelchair

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

TRANSACTION TYPE:

☐ Order ☐ Quote

Date _____ Account # _____ P.O. # _____

Company Name _____ Phone # _____

Address _____ City _____ State _____ Zip Code _____ Country _____

Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

FRAME SELECTIONS

FRAME TYPE

☐ SOLARA3G Solara 3G Tilt-in-Space wheelchair \$2,982.00 E1161

A seating system must be purchased. No seating system is provided. 300 lb standard weight capacity.

See www.invacare.com/Solara3G for additional information when configuring this chair.

CUSTOM PACKAGES

☐ HDPG Heavy Duty Package (400 lbs. Weight Capacity) \$100.00 K0108
Reinforced axles and crossbraces. Weight limit of 400 lbs.

SUSPENSION SELECTION

NEW SUSPENSION

☐ DROCK DeGage Dynamic Rocker Back \$1295.00 K0108

Not available with HDPG, TRRO, LOSHEAR, LOSHEARWO, RCANE, FLDBACK. See www.degage.us for more details.

☐ FLSF Frog Legs Suspension Forks \$379.00 E1015

Not available with HDPG, TRRO. See chart page 7.

☐ SPACK Suspension Package \$1449.00 E1015

Not available with HDPG, TRRO, LOSHEAR, LOSHEARWO, RCANE, FLDBACK. Includes DROCK and FLSF.

Transport Tie Down

TRRO includes four factory-installed transport brackets and wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNA WC Vol 1. Section 19 Frontal Impact Test requirements for wheelchairs with a 168 lbs. crash dummy, which corresponds to a person with a weight of 114 lbs. to 209 lbs.

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

TRANSPORT/VENT TRAY SELECTION

TRANSPORT TIE DOWN

☐ Clear

☐ TRRO Transport Ready Option \$250.00 K0108

Not available with W23, W24, HDPG, AVENT, BACKVENT, LOSHEAR, LOSHEARWO, P750, 2218, OFFMNT, BILAT, 110CP.

☐ TRBKTS Transport Brackets \$250.00 K0108

Not available with AVENT, BACKVENT, P750, 2218.

Now available as a service part kit for existing models.

VENT TRAY OPTION

☐ AVENT Ventilator Tray (each) \$295.00 E1029

Not available with HDPG.

Available with 16" wide frame and above and 18" Seat-to-Floor Height and above.

☐ ATBOX Ventilator Battery Box (each) \$97.00

☐ BACKVENT Laptop Vent Tray \$1,540.00 E1030

Not available with HDPG, LOSHEAR, P750.

Available with 15" wide frame and above.

FRAME SELECTIONS

TILTING OPTIONS

☐ FRLT Cable-Free Foot Release-Left Mount No Charge

☐ FRRT Cable-Free Foot Release-Right Mount No Charge

☐ THR Traditional Hand Release Standard

☐ TLL Tilt Lock Levers No Charge

☐ TLC Tilt Lock Stops (pair) \$30.00

Additional pair of tilt lock stops, one pair comes standard with chair.

SEAT WIDTH RANGE

☐ SWR1216 12" - 16" Seat Width Range No Charge

Not available with HDPG.

☐ SWR1418 14" - 18" Seat Width Range No Charge

☐ SWR1620 16" - 20" Seat Width Range No Charge

☐ SWR1822 18" - 22" Seat Width Range No Charge

☐ SWR2024 20" - 24" Seat Width Range No Charge

NOTE: All specifications and dimensions are approximate.

For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Solara™ 3G - Adult Tilt-in-Space Wheelchair

SEAT WIDTH

<input type="radio"/> W12	12" Wide	No Charge	
<input type="radio"/> W13	13" Wide	No Charge	
<input type="radio"/> W14	14" Wide	No Charge	
<input type="radio"/> W15	15" Wide	No Charge	
<input type="radio"/> W16	16" Wide	No Charge	
<input type="radio"/> W17	17" Wide	No Charge	
<input type="radio"/> W18	18" Wide	No Charge	
<input type="radio"/> W19	19" Wide	No Charge	
<input type="radio"/> W20	20" Wide	\$300.00	E2201
<input type="radio"/> W21	21" Wide	\$300.00	E2201 K0108
<input type="radio"/> W22	22" Wide	\$300.00	E2201 K0108
<input type="radio"/> W23	23" Wide	\$300.00	E2201 K0108
<input type="radio"/> W24	24" Wide	\$300.00	E2202 K0108

SEAT DEPTH RANGE

- ☐ SDR1216 12" – 16" Seat Depth Range No Charge
Tilt Assist feature not available on this frame.
- ☐ SDR1519 15" – 19" Seat Depth Range No Charge
- ☐ SDR1822 18" – 22" Seat Depth Range No Charge

SEAT DEPTH

<input type="radio"/> D12	12" Seat Depth	No Charge	
<input type="radio"/> D13	13" Seat Depth	No Charge	
<input type="radio"/> D14	14" Seat Depth	No Charge	
<input type="radio"/> D15	15" Seat Depth	No Charge	
<input type="radio"/> D16	16" Seat Depth	No Charge	
<input type="radio"/> D17	17" Seat Depth	No Charge	
<input type="radio"/> D18	18" Seat Depth	No Charge	
<input type="radio"/> D19	19" Seat Depth	No Charge	
<input type="radio"/> D20	20" Seat Depth	\$300.00	E2203
<input type="radio"/> D21	21" Seat Depth	\$300.00	E2203 K0108
<input type="radio"/> D22	22" Seat Depth	\$300.00	E2204 K0108

FRAME FINISH

FRAME COLORS

- ☐ 24P Wet Black No Charge
- ☐ 29P White Snow..... No Charge
- ☐ 60P Silver Metallic No Charge
- ☐ 115P Black Prism No Charge
- ☐ 105P Midnight Blue No Charge
- ☐ 116P Sunkissed Orange No Charge
- ☐ 117P Tangerine Red No Charge
- ☐ 30P Sunny Yellow No Charge
- ☐ 61P Electric Red No Charge
- ☐ 119P Bubblegum Pink No Charge
- ☐ 120P Grape Madness No Charge
- ☐ 62P Electric Blue No Charge
- ☐ 127P Grasshopper Green No Charge
- ☐ Touch Up Paint \$20.00
Touch up paint is not available in all colors.
- Touch Up Color

Discontinued frame colors are available through our specials department.
C-070109030544

CASTER, WHEEL AND TIRE SELECTIONS

SEAT-TO-FLOOR HEIGHT

- ☐ STF125 Seat-To-Floor 12 1/2" (10 1/2" w/DPSEAT) No Charge
- ☐ STF13 Seat-To-Floor 13" (11" w/DPSEAT) No Charge
- ☐ STF135 Seat-To-Floor 13 1/2" (11 1/2" w/DPSEAT) No Charge
- ☐ STF14 Seat-To-Floor 14" (12" w/DPSEAT) No Charge
- ☐ STF145 Seat-To-Floor 14 1/2" (12 1/2" w/DPSEAT) No Charge
- ☐ STF15 Seat-To-Floor 15" (13" w/DPSEAT) No Charge
- ☐ STF155 Seat-To-Floor 15 1/2" (13 1/2" w/DPSEAT) No Charge
- ☐ STF16 Seat-To-Floor 16" (14" w/DPSEAT) No Charge
- ☐ STF165 Seat-To-Floor 16 1/2" (14 1/2" w/DPSEAT) No Charge
- ☐ STF17 Seat-To-Floor 17" (15" w/DPSEAT) No Charge
- ☐ STF175 Seat-To-Floor 17 1/2" (15 1/2" w/DPSEAT) No Charge
- ☐ STF18 Seat-To-Floor 18" (16" w/DPSEAT) No Charge
- ☐ STF185 Seat-To-Floor 18 1/2" (16 1/2" w/DPSEAT) No Charge
- ☐ STF19 Seat-To-Floor 19" (17" w/DPSEAT) No Charge

CASTER OPTIONS (PAIR)

- ☐ 4SPLT 4" x 1" Composite Urethane Light Up No Charge
Not available with HDPG.
- ☐ 4ALUM2 4" x 1.4" Aluminum Soft Roll \$90.00
- ☐ 5SPLT 5" x 1 Composite Urethane Light Up No Charge
Not available with HDPG.
- ☐ 51DC 5" x 1" Urethane Casters No Charge
Not available with HDPG.
- ☐ 5ALUM2 5" x 1.4" Aluminum Soft Roll \$90.00
- ☐ 61DC 6" x 1" Urethane Caster No Charge
Not available with HDPG.
- ☐ 1259 6" x 2" Semi-Pneumatic Casters \$60.00
Not available with HDPG.
- ☐ 6ALUM2 6" x 1.4" Aluminum Soft Roll \$90.00
- ☐ 81DC 8" x 1" Urethane Casters No Charge
- ☐ 815DC 8" x 1.5" Urethane Caster No Charge

WHEEL OPTIONS (PAIR)

- ☐ 12COM 12" Composite Wheels No Charge
Not available with urethane tires.
- ☐ 16COM 16" Composite Wheels No Charge
Not available with urethane tires.
Only available with aluminum handrims.
- ☐ 18COM 18" Composite Wheels No Charge
Not available with urethane tires.
Only available with aluminum handrims.
- ☐ 20COM 20" Composite Wheels No Charge
- ☐ 22COM 22" Composite Wheels No Charge
- ☐ 24COM 24" Composite Wheels No Charge

AXLE TYPE (PAIR)

- ☐ A1577 Quick Release Axles No Charge
- ☐ PAXLE Permanent Axles Standard

TIRE OPTIONS (PAIR)

- ☐ PNEU Pneumatic Tires \$60.00
Not available with HDPG.
- ☐ PNEUFF Flat Free Inserts \$65.00 E2213
- ☐ URETH Urethane (Treaded) Tires No Charge

NOTE: All specifications and dimensions are approximate.

For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Solara™ 3G - Adult Tilt-in-Space Wheelchair

Front and Rear Seat-to-Floor Heights (measured from Floor to Top of Seat Rail)					
Rear Wheel	4" 4SPLT 4ALUM2	5" 51DC 5ALUM2 5SPLT	6" 61DC 6ALUM2	6" 1259	8" 81DC 815DC
12" Pneu PneuFF	12 1/2" 13" 13 1/2" 15 1/2" 16"	13 1/2" 15 1/2" 16" 16 1/2"	15 1/2" 16" 16 1/2" 17"	15 1/2" 16" 16 1/2" 17 1/2" 18"	16 1/2" 17" 17 1/2" 18"
12" Urethane	12 1/2" 14 1/2" 15" 15 1/2" 16"	14 1/2" 15" 15 1/2" 16" 16 1/2"	14 1/2" 15" 15 1/2" 16" 16 1/2"	14 1/2" 15" 15 1/2" 16" 16 1/2"	16 1/2" 17"
16"	12 1/2" 13" 13 1/2" 14" 14 1/2" 15"	13 1/2" 14" 14 1/2" 15"	14 1/2" 15" 15 1/2" 17"	14 1/2" 15" 17" 17 1/2" 18"	17" 17 1/2" 18" 18 1/2" 19"
18"	13 1/2" 14" 14 1/2" 15" 15 1/2" 16"	13 1/2" 14" 14 1/2" 15" 15 1/2" 16"	14 1/2" 15" 15 1/2" 16"	14 1/2" 15" 15 1/2" 16" 17" 17 1/2"	18" 18 1/2" 19"
20"	15" 15 1/2" 16"	15" 15 1/2" 16" 16 1/2"	15" 15 1/2" 16" 16 1/2" 17"	15" 15 1/2" 16" 16 1/2" 17" 17 1/2"	16 1/2" 17" 17 1/2" 18" 18 1/2"
22"	16"	16" 16 1/2"	16" 16 1/2" 17"	16" 16 1/2" 17" 17 1/2" 18" 18 1/2"	16 1/2" 17" 17 1/2" 18" 18 1/2"
24"	N/A	16 1/2"	16 1/2" 17"	16 1/2" 17" 17 1/2" 18" 18 1/2"	16 1/2" 17" 17 1/2" 18" 18 1/2" 19"
All Dimensions Are Approximate					

HANDRIM OPTIONS (PAIR)

- ☐ CHDRM Chrome Plated Handrims No Charge
☐ 1512 Plastic Coated Aluminum Handrims \$72.00
☐ AT-1225-8 Plastic Coated with 8 Vertical Projections \$240.00
Only available with 24COM7.
☐ AT-1224-8 Plastic Coated with 8 Oblique Projections \$240.00
Only available with 24COM7.
☐ AT1225-12 Plastic Coated with 12 Vertical Projections \$240.00
Only available with 24COM7.
☐ L8702 Omit Handrims Less \$15.00

WHEEL LOCKS

- ☐ FOOTBRK Foot Operated Push To Lock Wheel Lock \$50.00 K0108
Not available with AVENT.
☐ PULLLOCK Pull To Lock Wheel Locks (pair) No Charge
☐ PUSHLOCK Push To Lock Wheel Locks (pair) No Charge
☐ EXTLOCK Push To Lock with 6 Bolt-On Ext (pair) \$65.00 E0961
☐ FOOTHUB Foot Operated Hublock \$200.00 K0108
Not available with 16COM.
☐ ATHUB Cane Mount Hub Lock \$50.00 K0108
Not available with LOSHEAR or LOSHEARWO, 16COM.
☐ USEHUB Seat Rail Mount Hub Lock \$50.00 K0108
Only available with TAD of TAF.
Not available with 16COM.

ARM AND BACK SELECTIONS



SEAT STYLE

☐ Clear

FLSEAT or DPSEAT are solid seat pans and do not come standard with chair. Both pans grow in width and require no new hardware or parts when converting from flush to drop or vice versa.

- ☐ FLSEAT Flush Seat Pan \$95.00 E2231
☐ DPSEAT Drop Seat Pan \$95.00 E2231
Accommodates 2" lower seat to floor heights. Use a cushion not less than 3" thick.
☐ NOSEAT No Seat Pan \$0.00

Seat Pan Growth Range (1" increments)	
Chair Width	Seat Pan Width Range
W12	12"-14"
W13	12"-14"
W14	13"-16"
W15	13"-16"
W16	15"-18"
W17	15"-18"
W18	17"-20"
W19	17"-20"
W20	19"-22"
W21	19"-22"
W22	21"-24"
W23	21"-24"
W24	23"-24"

Arm Heights (measured from Top of Seat Rail to Top of Armpad)		
Arm Type	Low	High
TAD TAF	6 1/2" - 12 1/2" 1/2" increments	9" - 17" 1" increments
FXDESK FXFULL	9 1/2"	12"
ADJDESK ADJFULL	9 1/2" - 13 1/2" 1" increments	12" - 16" 1" increments
FADJDESK FADJFULL	9" - 13" 1" increments	12" - 16" 1" increments

ARM HEIGHT

☐ Clear

- ☐ SARMHT Low Arm Height No Charge
☐ HARMHT High Arm Height No Charge

CANTILEVER ARMS

☐ Clear

Not available with LOSHEAR or LOSHEARWO

- ☐ AT99N Cantilever Non Locking Bolt on Arms No Charge
☐ AT55D Cantilever Locking Desk Length Bolt on Arms \$195.00
☐ AT55F Cantilever Locking Full Length Bolt on Arms \$195.00

NOTE: All specifications and dimensions are approximate.

For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Solara™ 3G - Adult Tilt-in-Space Wheelchair

T-ARMS (PAIR)

☐ Clear

- ☐ TAD Adjustable Height, Desk Length \$225.00 E0973
- ☐ TAF Adjustable Height, Full Length \$225.00 E0973

TWO-POINT, FIXED HEIGHT ARMS (PAIR)

☐ Clear

Available on 15" seat depth and above.

- ☐ FXDESK Conventional, Desk Length No Charge
- ☐ FXFULL Conventional, Full Length No Charge

TWO-POINT, ADJUSTABLE HEIGHT ARMS (PAIR)

☐ Clear

Available on 15" seat depth and above.

- ☐ ADJDESK Conventional, Desk Length \$195.00 E0973
- ☐ ADJFULL Conventional, Full Length \$195.00 E0973
- ☐ FADJDESK Flip-Back, Desk Length \$245.00 E0973

Not available with LOSHEAR, LOSHEARWO or HDPG.

- ☐ FADJFULL Flip-Back, Full Length \$245.00 E0973

Not available with LOSHEAR, LOSHEARWO or HDPG.

BACK TYPE

- ☐ SCANE Standard Aluminum Back Cane Standard
- ☐ RCANE Steel-Reinforced Back Cane. No Charge
- ☐ LOSHEAR Low Shear Reclining Back with Upholstery \$810.00 K0108
- Not available with HDPG.
Comes with back upholstery and AASH.
Available on 14"-20" seat width.
- ☐ LOSHEARWO Low Shear Reclining Back without Upholstery \$790.00 K0108
- Not available with HDPG, headrest and pillow. Comes standard with AASH.
- ☐ RECEXT Reclining Headrest Pillow and Extension \$75.00 E0966
- ☐ RLL Recliner Lockout Levers No Charge
- Not available with SCANE or RCANE.
- ☐ RLC Recliner Lock Stops \$30.00

Back Cane Heights (measured from Top of Seat Rail)

Back Cane Type	Inches to 10 degree bend	Total Height
BH20	Straight Cane	20.5"
BH20B/BH24B	15 1/2"	24.5"
BH24	Straight Cane	24 1/2"
LOSHEAR	Straight Cane	26"
LOSHEAR + Headrest	Straight Cane	36"
ABH1720	16.2" - 21.2"	29.3" - 34.3"

BACK HEIGHT

☐ Clear

- ☐ BH20B Back Height 20" with bend No Charge
- Not available with P750.
- ☐ BH20 Strt Cane w Push Handle Back Ht 20" No Charge A9900
- Not available with P750.
- ☐ BH24B Back Height 24" with bend No Charge
- ☐ BH24 Back Height 24" No Charge
- ☐ ABH1720 Adj Height Back Canes (17" - 20") No Charge
- Not available with RCANE (because it consists of heavy duty construction), 1483, ATHUB, P750.

BACK CANE OPTIONS

☐ Clear

- ☐ FXBACK Non-Fold Down Back Canes Standard
- ☐ FLDBACK Fold-Down Back Canes No Charge
- Not available with LOSHEAR, DROCK.

POSITIONING STRAPS (EACH)

☐ Clear

- ☐ 1515 Push Button Buckle Seat Positioning Strap \$50.00 E0978
- ☐ 1321BK Hook and Loop Seat Positioning Strap \$50.00 E0978
- ☐ 1311BK Airline Buckle Seat Positioning Strap \$50.00 E0978
- ☐ 9515 Push Button Chest Positioning Strap \$100.00 E0960
- Not available with LOSHEAR or LOSHEARWO.
- ☐ 7321BK Hook and Loop Chest Positioning Strap \$100.00 E0960
- Not available with LOSHEAR or LOSHEARWO.
- ☐ 7311BK Airline Buckle Chest Positioning Strap \$100.00 E0960
- Not available with LOSHEAR or LOSHEARWO.

BODYPOINT® PADDED CALF STRAPS

☐ Clear

See www.bodypoint.com for complete details. Measure chair width to determine size.

- ☐ BPCALFS Padded Calf Straps (13" - 15") \$52.00 K0038
- ☐ BPCALFM Padded Calf Straps (15" - 18") \$56.00 K0038
- ☐ BPCALFL Padded Calf Straps (18" - 22") \$64.00 K0038

BODYPOINT® ANKLE HUGGERS®

☐ Clear

See www.bodypoint.com for complete details. Measure circumference above ankle to determine size.

- ☐ BPAHGS Side - Release Ankle Huggers (6.5" - 8") \$156.00 K0108
- ☐ BPAHGM Side - Release Ankle Huggers (7.5" - 9") \$160.00 K0108
- ☐ BPAHGL Side - Release Ankle Huggers (8.5" - 11.5") \$160.00 K0108
- ☐ BPAHGX Side - Release Ankle Huggers (11" - 13") \$164.00 K0108

BODYPOINT® STAYFLEX™

☐ Clear

See www.bodypoint.com for complete details. Measure shoulder width to determine size. BPSFM and BPSL are approximately 2" wider in center than BPSFNM and BPSFNL.

- ☐ BPSFM Anterior Trunk Support with Zip (13" - 16") \$174.00 E0960
- ☐ BPSFL Anterior Trunk Support with Zip (16" - 19") \$178.00 E0960
- ☐ BPSFNM Anterior Trunk Support Narrow with Zip (13" - 16") .. \$174.00 E0960
- ☐ BPSFNL Anterior Trunk Support Narrow with Zip (16" - 19") .. \$178.00 E0960

MISCELLANEOUS ACCESSORIES

BACK ACCESSORIES

- ☐ ADJRGD Adjustable Width Back Rigidizing Bar \$85.00 K0108
- ☐ 1483 Non Adjustable Stroller Handle \$85.00
- ☐ AASH Adj.Angle Stroller Handles (Outside Mount) No Charge
- ☐ AASHIN Adj.Angle Stroller Handles (Inside Mount) No Charge

MISCELLANEOUS ACCESSORIES


- ☐ 1209 Tire Pump \$35.00 A9270
- ☐ P750 Oxygen Holder \$128.00 E2208
- Not available with HDPG, LOSHEAR, LOSHEARWO, BH20, BH20B, ABH1720.
- ☐ 2218 IV Pole \$150.00 K0105
- Not available with LOSHEAR or LOSHEARWO.
- ☐ A1580 Rear Anti-Tippers No Charge
- ☐ BP02 Multi Pocket Backpack \$85.00 A9270

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

FRONT RIGGING OPTIONS

NEW FRONT RIGGING MOUNTING

- ☐ STDMNT STD Front Rigging Mount Standard
- ☐ OFFMNT Offset Front Rigging Mount No Charge
Not available with 12" - 16" seat depth range.



FRONT RIGGINGS (PAIR)

- ☐ 70LIFTNT 70° Lift Off Non-Taper Footrests \$95.00
- ☐ 70LIFTMFX 70° Lift Off MFX Footrests \$95.00
Must select AT5543.
- ☐ 90LIFTMFX 90° Lift Off MFX Footrests \$95.00
Must select AT5543.
- ☐ SWING1319 60° Swingaway Footrests No Charge
- ☐ SWING914 60° Swingaway Footrests MFX No Charge E0990
- ☐ 70NTAPER 70° Swingaway Footrests No Charge
- ☐ 70MFX 70° Swingaway Footrests MFX No Charge E0990
Must select AT5543.
- ☐ 70STAPER 70° Tapered Footrests No Charge E0990
Available on seat widths 16" and above.
- ☐ 80NTAPER 80° Swingaway Footrests No Charge
- ☐ 80MFX 80° Swingaway Footrests MFX No Charge
Must select AT5543.
- ☐ 90SWING 90° Swingaway Footrests No Charge
Must select AT5543.
- ☐ AT5044LIFT Lift Off Elevating Legrests \$255.00 E0990
- ☐ AT5044 Elevating Legrests \$255.00 E0990
- ☐ BILAT Bi-Lateral Contracture Footrests \$120.00
Not available with HDPG OR 12 1/2" STF Height.
- ☐ ATL4A Center Pivot Articulating Legrests \$320.00 K0053
- ☐ 110CP Contracture Platform \$280.00
Not available with HDPG.
Available on 16"-22" seat widths.

FOOTPLATE TYPE (PAIR)

- ☐ 1651 Flip-Up Composite Footplates No Charge
Not available with Bodypoint Ankle Huggers.
- ☐ A1562 Aluminum Footplates \$59.00
Not available with Bodypoint Ankle Huggers.
- ☐ A1564 Extra Large Aluminum Footplates \$59.00 K0041
Not available with Bodypoint Ankle Huggers.
- ☐ AT5543 Adjustable Angle Flip-Up Footplates \$165.00 K0040

ADJUSTABLE ANGLE FOOTPLATE RIGHT (EACH) (BI-LATERAL CONTRACTURE FOOTREST ONLY)

Must select with BILAT.

- ☐ AAF325R Adj. Angle Footplate Right 3 1/4" W \$83.00 K0040
- ☐ AAF425R Adj. Angle Footplate Right 4 1/4" W \$83.00 K0040
- ☐ AAF525R Adj. Angle Footplate Right 5 1/4" W \$83.00 K0040
- ☐ AAF625R Adj. Angle Footplate Right 6 1/4" W \$83.00 K0040
- ☐ AAF725R Adj. Angle Footplate Right 7 1/4" W \$83.00 K0040

ADJUSTABLE ANGLE FOOTPLATE LEFT (EACH) (BI-LATERAL CONTRACTURE FOOTREST ONLY)

Must select with BILAT.

- ☐ AAF325L Adj. Angle Footplate Left 3 1/4" W \$83.00 K0040
- ☐ AAF425L Adj. Angle Footplate Left 4 1/4" W \$83.00 K0040
- ☐ AAF525L Adj. Angle Footplate Left 5 1/4" W \$83.00 K0040
- ☐ AAF625L Adj. Angle Footplate Left 6 1/4" W \$83.00 K0040
- ☐ AAF725L Adj. Angle Footplate Left 7 1/4" W \$83.00 K0040

Knee-to-Heel Range (Measured from edge of seat pan to back of footplate)		
Front Rigging Type	Minimum	Maximum
Range adjusts in 1/2" increments.		
60°	13 1/2"	17 1/2"
60° MFX	10"	13 1/2"
*70°	13 1/2"	18 1/2"
*70° Tapered	13 1/4"	17.00"
80°	10.5"	16"
Bilateral Contracture STF125-STF16	5 1/2"	9"
Bilateral Contracture STF165-STF19	5 1/2"	13"
Range adjusts in 1" increments.		
70° MFX	6 1/2"	12 1/2"
80° MFX	4.5"	9"
90°	4 1/2"	7 1/2"
Infinite Range Adjustments		
*Elevating	14.00"	20.00"
Articulating	16.00"	24.00"
Contracture Platform STF125-STF16	7 1/4"	8 1/4"
Contracture Platform STF165-STF19	11 1/2"	16"
* ALPT shifts range by 2"		

FOOTREST/LEGREST ACCESSORIES (PAIR)

- ☐ 1337 Calf Strap \$24.00 K0038
- ☐ 9342 "H" Strap \$56.00 K0039
- ☐ FRMPRB ChairArmor Blue Neoprene Frame Protectors ... \$65.00 A9270
Only available with 70LIFTNT, SWING1319, 70NTAPER, 70STAPER
- ☐ FRMPRS ChairArmor Silver Neoprene Frame Protectors . \$65.00 A9270
Only available with 70LIFTNT, SWING1319, 70NTAPER, 70STAPER
- ☐ FRMPRR ChairArmor Red Neoprene Frame Protectors \$65.00 A9270
Only available with 70LIFTNT, SWING1319, 70NTAPER, 70STAPER
- ☐ FRMPRBK ChairArmor Black Neoprene Frame Protectors . \$65.00 A9270
Only available with 70LIFTNT, SWING1319, 70NTAPER, 70STAPER
- ☐ ALPT Longer Pivot And Slide Tubes \$84.00
Only available with 70LIFTNT, SWING1319, 70NTAPER, 70STAPER
- ☐ ALPT3 3" Extension For 90SWING Footrests \$31.00 K0108
- ☐ SCLAMP Clamp-On Sector Blocks..... \$50.00 K0108
Not available with 70LIFTNT, 70LIFTMFX, 90LIFTMFX, AT5044LIFT, BILAT, 110CP. Adjustable angle. See Front Rigging Mounts image.
- ☐ CLFPDLRG Large Calf Pad for Elevating Leg Rests.. No Charge
Only available on Elevating leg rests.

HEEL AND TOE LOOP (PAIR)

- ☐ A1600BK Heel Loops With Ankle Straps \$51.00 E0951
- ☐ HLOOP Heel Loops Without Ankle Straps \$51.00 E0951
- ☐ 1345 Toe Loops \$39.00 E0952

NOTE: All specifications and dimensions are approximate.

For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Solara™ 3G - Adult Tilt-in-Space Wheelchair

SEATING

SEATING CUSHION SELECTION

☐ Clear

<input type="radio"/> ITSG	Matrx Stabilite Gentle Contour Cushion	\$320.00	E2605 E2606
<input type="radio"/> ITFG	Matrx Flovair Gentle Contour Cushion	\$390.00	E2607 E2608
<input type="radio"/> ITFM	Matrx Flovair Max Contour Cushion	\$390.00	E2607 E2608
<input type="radio"/> MA-VI	Matrx Vi Cushion	\$425.00	E2607 E2608
<input type="radio"/> PS-V	Matrx PS Cushion	\$345.00	E2605 E2606
<input type="radio"/> FT	Matrx Flo-Tech Cushion	\$475.00	E2622 E2623

SEATING BACK SELECTION

☐ Clear

<input type="radio"/> MPB-V	Matrx PB Back	\$625.00	E2615
<input type="radio"/> MPBD	Matrx PB Deep Back	\$645.00	E2615
<input type="radio"/> PBE	Matrx Elite Back	\$625.00	E2620
<input type="radio"/> EDB	Matrx Elite Deep Back	\$760.00	E2620
<input type="radio"/> HDB	Matrx PB Heavy Duty Back	\$835.00	E2615 E2616
<input type="radio"/> DHDB	Matrx PB Deep Heavy Duty Back	\$935.00	E2615 E2616

HEADREST OPTIONS

☐ Clear

<input type="radio"/> CP175	Stealth Headrest 14" Multi Axis	\$225.00	E0955
<input type="radio"/> CP175HW	Removable and Adjustable Hardware	\$225.00	E1028
<input type="radio"/> CP275	Stealth Headrest 10" Multi Axis	\$225.00	E0955
<input type="radio"/> CP275HW	Removable and Adjustable Hardware	\$225.00	E1028
<input type="radio"/> CP180	Stealth Headrest 14" TWB	\$225.00	E0955
<input type="radio"/> CP180HW	Removable and Adjustable Hardware	\$225.00	E1028
<input type="radio"/> CP280	Stealth Headrest 10" TWB	\$225.00	E0955
<input type="radio"/> CP280HW	Removable and Adjustable Hardware	\$225.00	E1028

FROGLEGS SUSPENSION FORKS

Front & Rear-Seat-to-Floor Heights Measured from floor to top of seat rail					
Rear Wheel	4" 4SPLT	5" 5IDC 5SPLT 5ALUM2	6" 6IDC 6ALUM2	6" 1259	8" 8IDC 8I5DC
12"	13.50"	15.50"	15.50" 16.00"	16.50" 17.00" 17.50" 18.00"	N/A
16"	13.50" 14.00" 14.50" 15.00"	14.50" 15.00"	N/A	17.00" 17.50" 18.00"	18.50" 19.00"
18"	13.50" 14.00" 14.50" 15.00"	14.50" 15.00" 15.50"	15.50" 16.00"	18.00"	18.50" 19.00"
20"	15.00"	15.00" 15.50"	15.50" 16.00"	16.50" 17.00" 17.50"	N/A
22"	N/A	N/A	16.00"	16.50" 17.00" 17.50" 18.00"	18.50"
24"	N/A	N/A	N/A	16.50" 17.00" 17.50" 18.00"	18.50" 19.00"

Basic Chair Includes:

- Solara Tilt-in-Space Frame
- Seat-to-Floor Height Selection
- Tilt Option (FRLT, FRRT, THR)
- Casters (4SPLT, 5SPLT, 5IDC, 6IDC, 8IDC)
- Frame Width (12"-19")
- Wheels
- Axles (A1577, PAXLE)
- Urethane Treaded Tires
- Handrims (AHDRM, CHDRM)
- Center-of-Gravity Setting - Based on seat depth and seat depth range. See page 7 for details.
- Frame Finish
- Two-Point Conventional Fixed Arms (FXDESK or FXFULL)
- Rear Anti-Tippers
- Aluminum or Steel Reinforced Back Canes
- Wheel Locks (PULLLOCK, PUSHLOCK)
- Back Height Selection
- Front Riggings (SWING1319, SWING914, 70NTAPER, 70MFX, 70STAPER, 90SWING)
- Fold Down Back
- Flip-up Composite Footplates

NOTE: All specifications and dimensions are approximate.

For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Solara™ 3G - Adult Tilt-in-Space Wheelchair

Specifications	
Frame Style	Tilt-in-Space
Overall Width	Seat width + 8 1/4" (12" wheels) Seat width + 9 1/4" (18-24" wheels)
Seat Widths	12" - 24"
Seat Depths	12" - 22"
Back Height	17" - 20" , 20" , 24"
Adjustable Angle Back	0°, 5°, 10°, 15°, 20°, 25°, 30°
Tilt Range	-5° to 50°
Recliner Back	175°
Weight Limit	300 lbs./400 lbs. with Heavy-Duty Package
Weight	34 1/2 lbs*
Warranty	Lifetime on frame and crossbraces
*Base chair without Accessories	

Available Through Our Specials Departments

- 20" High Reclining Back SM-050503023551
 - 25" - 26" Wide SM-010312110036
 - 80° Footrest Hanger SM-200411121756
 - Front Rigging Mount: C-090731110749
- Standard on one side and offset on the opposite side

Tilt Assist Explained

The Solara® 3G wheelchair is shipped from the factory pre-configured with one of three "Tilt Assist" settings: *Front, Mid or Rear*, based solely on the seat depth range and seat depth ordered. This setting describes how the user's weight is distributed on the frame itself, i.e., does one front load the chair, rear load or is the weight evenly distributed. As compared to the 2G, the 3G requires up to 30% less force needed to tilt because of these settings. If factors other than seat depth or seat depth range come into play that affect the weight distribution of the chair, the provider has the ability to change the Tilt Assist setting in the field by simply adjusting 2 screws on each side of the frame.

Tilt Assist Factory Setting

	Seat Depth Range 15" - 19"	Seat Depth Range 18" - 22"
D15	Front Setting	
D16	Front Setting	
D17	Mid Setting	
D18	Rear Setting	Front Setting
D19	Rear Setting	Front Setting
D20		Mid Setting
D21		Rear Setting
D22		Rear Setting


The factory setting may need to be adjusted if a seating system is added. The factory setting does not account for conditions outside of the seat depth and seat depth range which may cause more weight to be distributed to the front or back of the chair.

When adjusting the Tilt Assist there are 4 bolts to change (2 on each side.) See the illustration to the left and chart above.

The **Rear** setting has the rocker tube telescoped entirely into the seat frame support tube where the adjustment bolt is in the bottom hole. Along the base frame the bolt is in the top hole. (e.g. accommodates a seat depth that is in the highest end of the seat depth range)

For the **Mid** setting the adjustment bolt is simply in the middle hole along the rocker tube which telescopes out and the middle hole along the base frame. (e.g. optimal for a seat depth in the middle of the seat depth range)

The **Front** setting has the bolt in the top hole of the rocker tube which is fully telescoped out of the seat frame support tube and in the bottom hole along the base frame. (e.g. accommodates a seat depth in the lowest end of the seat depth range)



Comments

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare Corporation

www.invacare.com

USA

One Invacare Way
Elyria, OH USA
44035-4190
800-333-6900

©2013 Invacare Corporation

All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.

