LONG TERM CARE BED ULTRA LOW & HIGH

★★★★★ 36" FIXED WIDTH

ULTRA L W BED



8 Function Hand Pendant



HI/LO, HEAD, FOOT, AUTO CONTOUR

EXTENSIVE TRAVEL RANGE: 7" - 30"

EXPANDABLE LENGTH: 80" - 88"

LOW POSITION OF ONLY 7"! **IDEAL SOLUTION FOR FALL PATIENTS**





Mattress Retainers included



Bumper Bar Protects Walls

• 8 Function hand pendant with universal symbols

- Auto contour function
- Quick Release boards
- Head and Foot Boards included
- Bumper bar protects walls from scratches and scuff marks
- Durable reinforced steel tube frame
- IV sockets on bed corners
- Corner mattress retainers included
- 8 3" CASTERS: 4 locking and 4 swivel
- HEAVY WEIGHT CAPACITY: 500 Lbs.



Quick Release Head & Foot Boards



Model #	Description	Width	Length	Low Range	High Range	Weight Cap
ULB7/30	Ultra Low	36"	80" - 88"	7"	30"	500