

Powerbase Wheelchair



P301 Gemini

Medicare Code: K0825 P301-2HD C
 K0824 P301-2HD S
 K0831 P3011-2SSE C
 K0830 P3011-2SSE S
Weight capacity: 450lbs
Weight capacity: 400lbs
 (w/optional seat lift)

Features

- Unique suspension system and high-density foam seat provide a smooth and comfortable ride
- The sophisticated controller coupled with 14" aluminum cast drive wheels offer excellent maneuverability.
- Comfortable 22" wide captain seat has swivel and folding capabilities with adjustable head and backrest
- Aluminum foot plate with flip-up height and angle adjustability provides various comfort positions, width adjustable armrests to provide various comfort positions.
- Programmable controller with a wide range of selections to suit the needs and demands of different users.
- Colors available : Red, Blue.

Options

Elevating legrests. Oxygen tank holder. Cane & crutch holder. Walker holder. Rear foldable basket. Canopy. Cover. Pan seat. Seat cushion. Cup holder. Side bag. Power seat elevator.

Features

- Real mid-wheel drive design for tight space maneuverability.
- Dual in-line motors for enhanced efficiency, torque, range and performance.
- Full front and rear suspension for easy obstacle climb and enhanced performance on varied terrain.
- 12" knobby tires and optimum weight distribution provide superior traction.
- New PG VR-2 or Dynamic Shark fully programmable controller.
- Off board charger and 4-post seat with semi-reclining backrest.
- NF-22 batteries for maximum range.
- Swing back joystick mount is standard.
- No rear tipping on footplate standing (up to 300-lbs).
- Colors available : Red, Blue.

Options

Elevating legrests. Oxygen tank holder. Cane & crutch holder. Walker holder. Foldable rear basket. Canopy. Cover. Pan seat. Seat cushion. Cup holder. Side bag. Power seat elevator.

Features

- Heavy-duty Power Wheelchair with a unique and prominent design.
- Four seat posts provide a more stable and durable structure for the frame.
- Extra-wide seat and extra-high back provide great comfort for bariatric users.
- Super high power in-line motors help riders negotiate 10° uphill ramps easily.
- Colors available : Red, Blue.

Options

Elevating leg rests. Oxygen tank holder. Cane & crutch holder. Walker holder. Foldable rear basket. Canopy. Cover. Pan seat. Seat cushion. Cup holder. Side bag.



P327 Vision Super

Medicare Code: K0825 P327-2HD C
 K0824 P327-2HD S
 K0831 P3274-2SSE C
 K0830 P3274-2SSE S
Weight capacity: 450lbs
Weight capacity: 400lbs
 (w/optional seat lift)



P710 Atlantis 1

Medicare Code: K0827 P710 2VHD C
 K0826 P710-2VHD S
Weight capacity: 600lbs



P301

Technical Specs

max.speed : 5mph
 range : up to 20mi.
 turning radius : 21"
 ground clearance : 4"
 motor : DC 24V, 240W
 controller : Dynamic Shark 60A / PG VR2 60A
 battery : 12V / 22NF(50AH) x 2
 charger: 5A off-board.
 gradient : 10°
 caster: 9" foam filled tires.
 drive wheel: 14" foam filled tires.
 brake : intelligent, regenerative, electromagnetic brakes

Weight (lb)

batteries	capacity
w/o 165	450 (400 w/ seat lift)
w/ 242	

Dimensions (in)

base	battery/ea	seat	fl w/batt.	capacity
22	20	22-24	20	
overall dimensions				
41	45	24		



P327

Technical Specs

max.speed : 5 mph
 range : up to 20mi.
 turning radius : 20.7"
 ground clearance : 2.2"
 motor : DC 24V, 300W x 2pcs
 controller : Dynamic Shark 60A
 battery : 12V/22NF x 2pcs
 charger: 5A off- board.
 gradient : 8°
 caster : F: 6" / R: 6" solid tires
 drive wheel : 12" pneumatic tires
 brake : intelligent, regenerative electromagnetic brakes

Weight (lb)

base	battery/ea	seat	fl w/batt.	capacity
121	38	46	243	450 (400 w/ seat lift)

Dimensions (in)

overall length	overall height	overall width	seat
35.4"	48	27.2	22
overall dimensions			
20	20.5	16.9-19.3	22.3-23.2



P710

Technical Specs

max.speed : 5mph
 range : up to 32mi.
 turning radius : 31.5"
 ground clearance : 3"
 motor : DC 24V, 400W
 controller : Dynamic Shark 60A / PG VR2 60A
 battery : 12V / Group 24 x 2
 charger: 6A off- board.
 gradient : 10°
 caster : 9" pneumatic tires
 drive wheel : 12" pneumatic tires
 brake : intelligent, regenerative electromagnetic brakes

Weight (lb)

batteries	capacity
w/o 218	600
w/ 304	

Dimensions (in)

base	battery/ea	seat	fl w/batt.	capacity
24	20	23-24		
overall dimensions				
25	44	50	32	