




Power Mobility Products

P326A VISION SPORT

Mid-Wheel Drive Power Base Wheelchair



P326A VISION SPORT

Mid-Wheel Drive Power Base Wheelchair

□ Technical Specs :

max.speed : 5mph
 range : up to 18 mi.
 turning radius : 20"
 ground clearance : 2"
 motor : DC 24V, 160W
 controller : Dynamic Shark 50A
 battery : 12V / U1 x 2 pcs
 charger : 4A off-board
 gradient : 6°
 caster : F: 6" / R: 6" solid tire
 drive wheel : 10" solid tires
 brake : intelligent, regenerative
 electromagnetic brakes

□ Features

- . Full suspension.
- . No rear tipping on footplate standing (up to 250-lbs).
- . 300-lb weight capacity.
- . Dual in-line motors for enhanced efficiency, torque, range and performance.
- . Swing back joystick mount is standard.
- . Height and width adjustable arms.
- . Semi-reclining back.
- . 20"captain seat.


□ Options:

Elevating leg rests. Oxygen tank holder.
 Cane & crutch holder. Walker holder.
 Foldable rear basket. Canopy.

□ Weight (lb)

base	battery/ea	seat	tl w/ball.	capacity
90	26	42	184	300

□ Dimensions : (in)



Color fender accessories



Medicare Code:
 P326-2S S (K0822)
 P326-2S C (K0823)

