




# Power Mobility Products

P327 VISION SUPER Bariatric





# P327

## VISION SUPER Bariatric

The P327 mid-wheel drive delivers outstanding maneuverability in tight, compact spaces and a smooth and stable ride both indoors and outdoors even over rough and uneven terrain. Its high performance suspension, together with the advantages of in-line motor technology for enhanced efficiency, torque, range and performance make the P327 the best in its class.

### Technical Specs :

max.speed : 5 mph  
 range : up to 20 mi  
 turning radius : 20.7"  
 ground clearance : 2.2"  
 motor : DC 24V,300W x 2pcs  
 controller : Dynamic Shark 60A  
 battery : 12V/22NF x 2pcs  
 charger : 5A off-board  
 gradient : 8°  
 caster : F: 6" / R: 6" solid tires  
 drive wheel : 12" pneumatic tires  
 brake : intelligent, regenerative  
 electromagnetic brakes


### Options :

Elevating legrests. Oxygen tank holder.  
 Gel cushion.

### Weight : (lb)

base	battery/ea	seat	tl w/batt.	capacity
121	38	46	243	450

### Dimensions : (in)



35.4" | 48 | 27.2 | 22



20 | 20.5 | 16.9-19.3 | 22-23.2



Medicare Code:  
 K0825 P327-2HD C  
 K0824 P327-2HD S  
 Seat Elevator:  
 K0831 P3274-2SSE C  
 K0830 P3274-2SSE S

### Features :

- . Real mid-wheel drive design for tight space maneuverability
- . Dual in-line motors for enhanced efficiency, torque, range and performance
- . Full front and rear suspension for easy obstacle climb and enhanced performance on varied terrain
- . 12" knobby tires and optimum weight distribution provide superior traction
- . Dynamic Shark fully programmable controller
- . 4-post seat with semi-reclining backrest
- . NF-22 batteries for maximum range
- . Off board charger
- . Swing back joystick mount is standard
- . No rear tipping on footplate standing (up to 300-lbs)