

Midi Electric Scooter

S141 Pioneer 4

Technical Specs

max.speed : 5mph
range : up to 18mi.
turning radius : 52"
ground clearance : 4"
motor : DC 24V, 250W
controller : Dynamic Rhino 70A
battery : 12V / U1(34AH) x 2
charger : 5A off-board
gradient : 8°
wheels : 10" pneumatic
brake : intelligent, regenerative electromagnetic brakes

Features

- Tiller adjusts to fit various arm lengths and collapses for easy storage.
- Sliding captain seat.
- Flip-up armrests with height, width & angle adjustment.
- Headlight, tail lights, turn signals & hazard lights.
- Off-board charger.
- Limited lifetime warranty on frame.
- Colors available : Red, Blue.

Options


Oxygen tank holder. Cane & crutch holder.
Walker holder. Canopy. Foldable rear basket.

Weight (lb)



batteries	capacity
w/o 137	400
w/ 189	

Dimensions (in)



Medicare Code

K0899

