



P101

Technical Specs

max.speed : 4mph
 range : 18mi.
 turning radius : 23"
 ground clearance : 4.8"
 motor : DC 24V, 180W
 controller : PG VSI
 battery : 12V / U1(34AH) x 2pcs
 charger : 5A off-board
 gradient : 10°
 caster : 8" foam filled tires
 drive wheel : 12-1/2" foam filled tires
 brake : intelligent, regenerative, electromagnetic brakes

Weight (lb)

batteries	capacity
w/o 86	300
w/ 150	

Dimensions (in)

seat		arms	
18	16	20	8-12
16	16	20	8-12
20	18	20	8-12

back	footeasts	
15-19	16	21
15-19	16	21
15-19	16	21

overall dimensions			
43	36	24	13
43	36	22	13
45	36	26	14



P181

Technical Specs

max.speed : 4mph
 range : up to 25mi.
 turning radius : 25"
 ground clearance : 4.8"
 motor : DC 24V, 240W
 controller : PG VSI 70A
 battery : 12V / 22NF(50AH) x 2pcs
 charger : 5A off-board
 gradient : 8°
 caster : 8" foam filled tires
 drive wheel : 12-1/2" foam filled tires
 brake : intelligent, regenerative, electromagnetic brakes

Weight (lb)

batteries	capacity
w/o 104	450
w/ 182	

Dimensions (in)

seat		arms	
22	18	20	8-12

back	footeasts	
15-19	16	21

overall dimensions			
45	36	28	15



P182

Technical Specs

max.speed : 4mph
 range : up to 25mi.
 turning radius : 26"
 ground clearance : 4.8"
 motor : DC 24V, 400W
 controller : PG VSI 70A
 battery : 12V / 22NF(50AH) x 2pcs
 charger : 5A off-board
 gradient : 10°
 caster : 9" foam filled tires
 drive wheel : 12-1/2" foam filled tires
 brake : intelligent, regenerative, electromagnetic brakes

Weight (lb)

batteries	capacity
w/o 108	600
w/ 186	

Dimensions (in)

seat		arms	
24	20	21	8-12

back	footeasts	
15-19	16	21

overall dimensions			
46	36	30	15



P183

Technical Specs

max.speed : 4mph
 range : up to 25mi.
 turning radius : 27"
 ground clearance : 4.8"
 motor : DC 24V, 400W
 controller : PG VSI 70A
 battery : 12V / 22NF(50AH) x 2pcs
 charger : 5A off-board
 gradient : 10°
 caster : 9" foam filled tires
 drive wheel : 12-1/2" foam filled tires
 brake : intelligent, regenerative, electromagnetic brake

Weight (lb)

batteries	capacity
w/o 112	700
w/ 190	

Dimensions (in)

seat		arms	
26	20	21	8-12


back	footeasts	
15-19	16	21

overall dimensions			
47.5	36	32	15

Folding Power Wheelchair

Folding Power Wheelchair

Travel-Ease Commuter Series



Features

- Folding battery tray allows for improved portability.
- Programmable controller with a wide range of selections to suit the needs and demands of different users.
- Auxiliary brakes offer extra safety.
- Adjustable height flip-up armrests, swing away footrests.
- Limited lifetime warranty on frame.
- Colors available : Silvertone.

Options

Elevating legrests. Oxygen tank holder.
Gel cushion. Back cushion.



P181

Medicare Code: K0824
weight Capacity: 450lbs
Seat Width: 22"



P101

Medicare Code: K0822
weight Capacity: 300lbs
Seat Width: 16", 18" or 20"

Features

- Easily foldable battery bracket allows for improved portability.
- Lightweight and foldable for easy handling, storage and transport.
- Programmable controller with a wide range of selections to suit the needs and demands of different users.
- Rugged heavy-duty casters.
- Adjustable height flip-up armrests, swing away footrests.
- 23" turning radius allows the chair to be used in very confined spaces.
- Limited lifetime warranty on frame.
- Colors available : Silvertone.

Options

Elevating legrests. Oxygen tank holder.
Gel cushion. Back cushion. Swing away controller bracket.

Features

- Folding battery tray allows improved portability.
- Programmable controller with a wide range of selections to suit the needs and demands of different users.
- Auxiliary brakes offer extra safety.
- Adjustable height flip-up armrests, swing away footrests.
- Limited lifetime warranty on frame.
- Colors available : Silvertone.
- Swing away controller bracket.

Options

Elevating legrests. Oxygen tank holder.
Gel cushion. Back cushion.



P182

Medicare Code: K0826
weight Capacity: 600lbs
Seat Width: 24"

Features

- Folding battery tray allows for improved portability.
- Programmable controller with a wide range of selections to suit the needs and demands of different users.
- Auxiliary brakes offer extra safety.
- Adjustable height flip-up armrests, swing away footrests.
- Limited lifetime warranty on frame.
- Colors available : Silvertone.
- Swing away controller bracket.

Options



P183

Medicare Code: K0828
weight Capacity: 700lbs



MERITS HEALTH PRODUCTS, INC.
730 N.E. 19th Place Cape Coral, FL 33909
TEL: 1-800-963-7487, 1-239-772-0579
SALES FAX: 1-239-574-2661
MAIN FAX: 1-239-772-3252