




KALMIA™

PERFECT HEIGHT

BED SYSTEM

EASY TO REST, EASY TO RISE

Meet the Kalmia Perfect Height Bed System. Whether you're short, tall, or in between, this bed features an automated control – you can enter and exit at the optimum comfort level. It's got all the adjustable features you need: head and foot elevation, zero gravity to reduce muscle pain and sleep apnea, a massage feature, and a climate controlled mattress. It comes in multiple sizes and looks like any premium bed you'd be happy to show in your home. **It's just easy!**



WE'VE MARRIED CUTTING EDGE ELECTRONICS WITH SLEEK HOME DÉCOR FOR COMFORT, CARE AND CONVENIENCE.



No jumping up or bending down. This bed meets you where you stand. Here's what to expect from the Kalmia Perfect Height Bed System.™

Hi-Low Technology – the bed rises and lowers a full 10" to make getting in and out easy

Adjustable from head to toe – custom-fit for your body

Therapeutic features, fashionable design – discreet electronics

Massage features – fits with 110 or 220 volt motors

KALMIA THERAPEUTIC MATTRESS— THE NEW TECHNOLOGY OF SLEEP

Better circulation – helps reduce swelling in legs

Customized support – relieves muscle and joint pain

Elevation control – decreases acid reflux / snoring

Bio-sensitive mattress – lowers risk of bed sores

Micro-climate ticking – breathable mattress stabilizes temperature

Side guard – reduces risk of falling

Sturdy design – 800 lb. capacity / 400 lb. capacity for Twin XL

White glove delivery and assembly available – we do it all

	Size	Load Capacity	Total Weight
Twin	36 x 80 inches	400 lbs	350 lbs
Full XL	54 x 80 inches	500 lbs	460 lbs
Queen	60 x 80 inches	500 lbs	470 lbs
Split King (2 Twins)	72 x 80 inches	400 lbs Each Side	700 lbs

Motorized Height Range (Deck to Floor): 10-20 inches



Includes Easy Control Remote

- DC-powered with emergency power-down
- Compatible with 110 or 220 volts
- Control up or down for easy in and out of bed
- No programming required
- Convenient pin-light flashlight for finding things in the dark or navigating a dark room



Sleep 37.5™

Bedding that maintains a microclimate of personal comfort while you sleep.

Bedding made with 37.5 technology actively removes excess heat and moisture vapor on-demand to maintain optimum comfort zone, allowing you to enter deeper, longer periods of sleep. Every Kalmia mattress features this innovation.

The Kalmia Therapeutic Sleep System by Parks Health Products

Our products are informed by more than 100 years of experience. From the latest breakthroughs in foam and fiber to integrated assembly and therapeutic systems, we deliver the entire solution – Parks Health is the new technology of sleep.



KALMIA™
THERAPEUTIC SLEEP SYSTEM

Contact your dealer today.