

Jazzy®

ELITE ES-1

The Jazzy® Elite ES-1 offers a combination of great style, comfort and performance for your active lifestyle. The in-line, front-wheel drive design provides tight turns around corners for optimal indoor maneuverability in tight spaces. Large front wheels offer extra absorption over obstacles and superb climbing capabilities combine with front anti-tip wheels for transitioning over grade changes easily.

IN-LINE FRONT-WHEEL DRIVE TECHNOLOGY



PERFORMANCE-MATCHED ELECTRONICS



COMFORT HIGH-BACK SEAT WITH HEADREST

AMERICA'S FAVORITE POWER CHAIR

Learn more online.



gojazzy.com

Specifications

Weight Capacity	300 lbs.
Maximum Speed	Up to 4 mph
Ground Clearance	0.75"
Turning Radius	23.5"
Overall Length	34.75"
Overall Width	22.75"
Drive Wheels	9" solid
Front Wheels	5" solid anti-tips
Rear Wheels	6" solid casters
Suspension	Limited
Drivetrain	Two-motor, in-line, front-wheel drive
Braking System	Regenerative and electro-mechanical
Standard Electronics	40 amp, PG GC 3 controller
Battery Charger	Off-board, 3 amp
Per-Charge Range	Up to 12 miles
Battery Requirements	(2) 12 volt, deep cycle, U-1
Battery Weight	25 lbs. each
Base Weight	69 lbs.
Standard Seat Weight	43 lbs.
Standard Seat	18" x 18-20"
Max. Seat Dimensions	20" x 20"
Warranty	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

Colors



Candy
Apple Red



Viper
Blue

Features

- Color-through shroud
- Black vinyl, depth adjustable high-back seat
- Flip-up foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

Options

- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Lap belt (50", 70", 80", 90")
- Swing-away footrest
- Elevating leg rests
- Heel loops
- Military patches
- Swing-away joystick
- Height adjustable footrest bracket
- XLR USB charger

Recommended Lifts

- Backpacker® Plus
- Backpacker® AVP
- Outlander Full Platform
- Zeus
- Commander



Use the Pride Lift Finder to find a compatible lift for your mobility product:

pridemobility.com/Lift-Finder

Models

Jazzy Elite ES-11S-C | Code: K0816