




Protect your investment with the Pride® Power Lift Recliner Care Kit



LC-105 shown in Sandal Microsuede fabric

Essential COLLECTION

The Essential Collection from Pride® Power Lift Recliners provides exceptional value and outstanding comfort. Several styles and fabric options will match any décor while delivering rest and worry-free operation.

AMERICA'S #1
POWER LIFT RECLINER

POWER
RECLINERS

1-800-800-8586 (U.S.) 1-888-570-1113 (Canada)
pridemobility.com prideliftchairs.com

Pride
Mobility Products Corp.
— LIVE YOUR BEST —



LC-105 shown in
Sandal Microsuede fabric

LC-105 FABRIC SELECTION



MICRO-SUEDE



COCOA



SANDAL



SKY

LEXIS URETHANE®



NEW CHESTNUT



Look for the Quick Ship
Indicator for the shortest
lead times available.



3-POSITION
FULL RECLINE

The 3-position full recline is available
on most Pride® Power Lift Recliner
models and puts you in a position ideal
for napping or relaxing.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.
The information contained herein is correct at the time of publication; we reserve the right to
alter specifications without prior notice.

SPECIFICATIONS

	Medium
	LC-105
	3-Position Chaise Lounger 5'2" - 5'10"
Weight Capacity	325 lbs.
Seat-to-Floor	18"
Seat Depth	19"
Seat Width	20"
Top of Back to Seat	27.5"
Back Style	Sewn Pillow Back
Fully Padded Chaise	Yes
Heat & Massage Option	No
Total Comfort Seating	No
Distance From Wall	18"
Overall Chair Width	32.75"
Chair Weight	96 lbs.
Warranty	1-Year Limited

BENEFITS

- Wooden frame composed of engineered furniture grade laminate/hardwood for added strength and durability
- Quiet & smooth lift system
- Heavy-duty lift actuator & scissor mechanism
- Standard head & arm covers
- Fire retardant (California specifications; bulletin 117-2013)
- Standard, large dual pockets
- Removable back
- 325 lbs. weight capacity
- Integrated battery backup in case of power failure¹
- Transformer features self-diagnostic electronics with 9' cord
- Lifetime warranty on steel lift frame
- 1-year limited warranty

OPTIONS*

- Residential cleaning kit
- Protective pad
- Non-skid leg levelers

Additional lead time may
apply for some options.

*Additional charges apply for all options.

1. Batteries not included.