



LC-525i shown in Aria® Cool Grey fabric

MADE IN USA  
OF FOREIGN & DOMESTIC COMPONENTS

USB  
Charger  
Port



# INFINITY COLLECTION

The Infinity Collection from Pride® Power Lift Recliners features true infinite-position seating and pocketed coils designed to deliver comfort, style and performance. Available in four different fabric selections, the Infinity Collection models deliver maximum relaxation and will match any décor.

AMERICA'S #1  
POWER LIFT RECLINER

POWER  
RECLINERS

Protect your investment with the Pride® Lift Chair Care Kit



US: 1-800-800-8586 • CAN: 1-888-570-1113  
pridemobility.com |

**Pride**  
Mobility Products Corp.  
— LIVE YOUR BEST —

# INFINITY COLLECTION

AMERICA'S #1  
POWER LIFT RECLINER

POWER  
LIFT  
RECLINERS



LC-525i shown in Ultraleather™ Pecan fabric



|                            | Petite Wide   | Small   | Medium  | Large   |
|----------------------------|---|---|---|---|
| <b>SPECIFICATIONS</b>      | <b>LC-525iPW<sup>2</sup></b>  | <b>LC-525iS<sup>2</sup></b>   | <b>LC-525iM<sup>2</sup></b>   | <b>LC-525iL<sup>2</sup></b>   |
|                            | True Infinite-Position<br>Zero Gravity<br>Trendelenburg<br>Chaise Lounger<br>5'2" & Below | True Infinite-Position<br>Zero Gravity<br>Trendelenburg<br>Chaise Lounger<br>5'3" & Below | True Infinite-Position<br>Zero Gravity<br>Trendelenburg<br>Chaise Lounger<br>5'4"- 6'0" | True Infinite-Position<br>Zero Gravity<br>Trendelenburg<br>Chaise Lounger<br>5'9"- 6'2" |
| <b>Weight Capacity</b>     | 375 lbs.  | 375 lbs.  | 375 lbs.  | 375 lbs.  |
| <b>Seat-to-Floor</b>       | 19"   | 19"   | 20"   | 21"   |
| <b>Seat Depth</b>          | 19"   | 19"   | 20"   | 22"   |
| <b>Seat Width</b>          | 24"   | 19"   | 20"   | 22"   |
| <b>Top of Back to Seat</b> | 28"   | 26"   | 29"   | 30"   |
| <b>Back Style</b>          | Blown Fiber<br>Split Back   | Blown Fiber<br>Split Back   | Blown Fiber<br>Split Back   | Blown Fiber<br>Split Back   |
| <b>Fully Padded Chaise</b> | Yes   | Yes   | Yes   | Yes   |
| <b>Heat and Massage</b>    | Optional  | Optional  | Optional  | Optional  |
| <b>Seating Type</b>        | Pocketed Coil<br>Springs  | Pocketed Coil<br>Springs  | Pocketed Coil<br>Springs  | Pocketed Coil<br>Springs  |
| <b>Distance From Wall</b>  | 19"   | 20"   | 20"   | 21"   |
| <b>Overall Chair Width</b> | 38"   | 34"   | 34.5"   | 37"   |
| <b>Chair Weight</b>        | 130 lbs.  | 130 lbs.  | 130 lbs.  | 130 lbs.  |
| <b>Warranty</b>            | 7-Year Prorated*  | 7-Year Prorated*  | 7-Year Prorated*  | 7-Year Prorated*  |

## INFINITY FABRIC SELECTION

QS

### DURASOFT



### LEXIS STA-KLEEN® sta\*kleen. PERFORMANCE FABRIC



### ULTRALEATHER™ Uf



### ARIA<sup>3</sup> CRYPTON. super fabric



Look for the Quick Ship Indicator for the shortest lead times available.

## BENEFITS

- Premium furniture grade laminate/hardwood frames for added strength and durability
- Pocketed coils deliver individualized support and relaxation
- Industry-first hand control with USB charger port **(not available with heat and massage)**
- Standard lumbar pillow
- Standard head and arm covers
- Large, dual pockets
- Easy to remove back
- Integrated battery backup in case of power failure<sup>1</sup>
- Transformer features self-diagnostic electronics with 9' cord

## OPTIONS\*

- Deluxe heat and massage
- Residential cleaning kit
- Head pillow
- Protective pad
- Lift chair overlay

**Additional lead time may apply for some options.**

\*Additional charges apply for all options.

1. Batteries not included.
2. Available with deluxe heat and massage only. Standard heat and massage not available.
3. Crypton<sup>®</sup> fabric offers a five-year warranty on the moisture barrier. Does not cover punctures in the fabric.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.  
\*Lifetime warranty on steel lift frame.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

US: 1-800-800-8586 • CAN: 1-888-570-1113  
pridemobility.com |

**Pride**  
Mobility Products Corp.  
— LIVE YOUR BEST —