

AN INCREDIBLE CHAIR! I CAN FULLY RECLINE AND IT IS SO RELAXING. IT REALLY IS THE MOST COMFORTABLE RECLINER IN THE WORLD!

- GEORGE R.



SCOTT MEUSER CHAIRMAN & CEO PRIDE MOBILITY®

With over three decades in the power lift recliner industry, we know that our business is about you, our customer. We're humbled when you choose to bring one of our power lift recliners into your home, knowing that it must fit your lifestyle and decor.

VivaLift!® Is the culmination of your direct feedback. You wished for luxurious comfort, exceptional quality, state-of-the-art features, and elegant designs. Over a two-year period, we listened and refined, then listened some more. The result is the VivaLift!® line of power lift recliners designed for you, your home and your lifestyle.

We understand that bringing a power lift recliner into your home isn't just a purchase, but an investment – toward your comfort, home décor, and lifestyle. We deeply value that you, our customer, have enabled us to meet those vital goals through VivaLift!®.









2. The upward lift offers convenient stand assist for getting in and out of the recliner.

3. The power lumbar enables you to personalize your comfort.

4. The footrest extension allows you to relax and stretch out completely.













5. The backrest moves up and down independently to help you find the perfect position.




6. Infinite position recline enables the backrest and footrest to move independently from one another so that a variety of positions can be reached.




THE MOST COMFORTABLE RECLINER IN THE WORLD!



STANDARD FEATURES







ULTRA

Exclusive VivaLift!® Tilt reclines the lift chair into a position that provides the feeling of weightlessness.

Advanced massage system with three intensity settings provides gentle pressure to your lower back.

New easy-to-use toggle remote features a memory button, tilt button and USB charger.

Convenient hideaway wireless phone charger and cup holder are located in the armrests.

True-infinite positioning enables the backrest and footrest to move independently from one another.

Lithium battery backup ensures operation in the case of a power outage.

STANDARD FABRIC: CAPRICCIO

HEADREST MEMORY

HEADREST LUMBAR

UP / DOWN BACK

LEG REST TILT







FDA Class II Medical Device¹

RADIANCE

Heat pads in the shoulders, lumbar and seat can be adjusted independently to suit your desired level of heat.

New easy-to-use toggle remote features 2 memory buttons, heat button and USB charger.

Convenient hideaway USB charging port and cup holder are located in the armrests.

True-infinite positioning enables the backrest and footrest to move independently from one another.

Lithium battery backup ensures operation in the case of a power outage.

HEAT HEAT

Heat pads are located

in the shoulder, lumbar

STANDARD FABRIC: CANYON





ATLAS PLUS 2

The power headrest lifts your head, neck and shoulders for maximum support.

Easy-to-use toggle remote features 2 memory buttons and USB charger.

Convenient hideaway USB charging port and cup holder are located in the armrests.

True-infinite positioning enables the backrest and footrest to move independently from one another.

Lithium battery backup ensures operation in the case of a power outage.

STANDARD FABRIC: BADLANDS



TRANQUIL 2

The power headrest lifts your head, neck and shoulders for maximum support.

The power lumbar enables you to personalize your comfort.

True-infinite positioning enables the backrest and footrest to move independently from one another.

Easy-to-use toggle remote features 2 memory buttons and USB charger.

Lithium battery backup ensures operation in the case of a power outage.





ELEGANCE 2

The power headrest lifts your head, neck and shoulders for maximum support.

The power lumbar enables you to personalize your comfort.

True-infinite positioning enables the backrest and footrest to move independently from one another.

Easy-to-use toggle remote features 2 memory buttons and USB charger.

Lithium battery backup ensures operation in the case of a power outage.

STANDARD FABRIC: BADLANDS

FDA Class II Medical Device¹



METRO 2

The power headrest lifts your head, neck and shoulders for maximum support.

The power lumbar enables you to personalize your comfort.

True-infinite positioning enables the backrest and footrest to move independently from one another.

Easy-to-use toggle remote features 2 memory buttons and USB charger.

Lithium battery backup ensures operation in the case of a power outage.

nd ner.

STANDARD FABRIC: SAVILLE





COLLECTION	ULTRA ¹	ULTRA ¹	ULTRA ¹	RADIANCE ^{1,2,4}	RADIANCE ^{1,2,4}	RADIANCE ^{1,2,4}	RADIANCE ^{1,2,4}	ATLAS PLUS 21,4
MODEL	PLR-4955S	PLR-4955M	PLR-4955L	PLR-3955PW	PLR-3955S	PLR-3955M	PLR-3955LT	PLR-2985M
POSITION ⁴	True-Infinite	True-Infinite	True-Infinite	True-Infinite	True-Infinite	True-Infinite	True-Infinite	True-Infinite
HEIGHT RANGE	5'4" and below	5′5″ - 6′	6' and above	5′4″ and below	5'4" and below	5′4″ - 6′	6' and above	5′4″ - 5′9″
SIZE	Small	Medium	Large	Petite Wide	Small	Medium	Large/Tall	Medium
WEIGHT CAPACITY	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.
SEAT-TO-FLOOR	18″	19.5″	21.5″	19.5″	18″	19.5″	21.5″	19.5"
SEAT DEPTH	18″	21"	22"	18″	18"	21"	22"	21"
SEAT WIDTH	20″	22.5″	23″	24"	20″	22.5″	23.5″	22"
TOP OF BACK TO SEAT	25″	28″	30″	23″	25″	28″	31″	26"
FULLY PADDED CHAISE	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
FOOTREST EXTENSION	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard
HEAT	Standard	Standard	Standard	Standard	Standard	Standard	Standard	NA
MASSAGE	Standard	Standard	Standard	NA	NA	NA	NA	NA
SEATING TYPE	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort
CHAIR HEIGHT	41"	43"	47"	40"	41"	43"	49"	42"
DISTANCE FROM WALL	15″	15″	21.5″	14"	15"	15″	21.5″	16"
OVERALL WIDTH	37"	39"	39"	41"	37"	39"	40"	39"
CHAIR WEIGHT	158 lbs.	158 lbs.	176 lbs.	161 lbs.	156 lbs.	158 lbs.	178 lbs.	164 lbs.
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated

^{*}Not all units achieve all positions. Please refer to specific product.

Note: Measurements may vary \pm 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

COLLECTION	TRANQUIL 21	TRANQUIL 21	TRANQUIL 21	TRANQUIL 21	ELEGANCE 21	ELEGANCE 21	METRO 21	URBANA 21
MODEL	PLR-935PW	PLR-935S	PLR-935M	PLR-935LT	PLR-975M	PLR-975L	PLR-925M	PLR-965M
POSITION ⁴	True-Infinite							
HEIGHT RANGE	5′3″ and below	5′3″ and below	5′4″ - 5′9″	5'10" and above	5′4″-5′9″	5′10″ and above	5′3″ and below	5'3" and below
SIZE	Petite Wide	Small	Medium	Large/Tall	Medium	Large	Medium	Medium
WEIGHT CAPACITY	400 lbs							
SEAT-TO-FLOOR	18"	18"	19.5"	21.5"	20"	20.5"	19"	20"
SEAT DEPTH	18.5"	18.5"	20.5"	22.5"	21"	22"	21"	20"
SEAT WIDTH	24"	20"	22"	23.5"	21"	22"	21"	21"
TOP OF BACK TO SEAT	24.5"	25"	27"	30"	27"	29"	26.5"	25"
FULLY PADDED CHAISE	Yes							
FOOTREST EXTENSION	Standard							
HEAT	NA							
MASSAGE	NA							
SEATING TYPE	Total Comfort							
CHAIR HEIGHT	38"	40"	42"	47"	42"	43"	41"	42"
DISTANCE FROM WALL	12.5"	13"	16"	16"	14"	16"	14"	14"
OVERALL WIDTH	36"	33.75"	37.5"	38"	35.5"	37.5"	33"	33"
CHAIR WEIGHT	143 lbs.	137 lbs.	145 lbs.	155 lbs.	152 lbs.	157 lbs.	142 lbs.	145 lbs.
WARRANTY	7-Yr. Pro-rated							

^{1.} Vivalift!® Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

^{2.} The wireless charger is only compatible with smart phones that have wireless charging capabilities. Additionally, no other smart devices, (i.e. smart watches) should be used with the wireless charger.

^{3.} Height ranges should be used as a guide only due to varying factors of each individual user.

^{4.} Unit has three different heating pads in the shoulders, lumbar and seat.



AUTHORIZED PRIDE® PROVIDER



pridemobility.com





