

**VivaLIFT!**<sup>®</sup>  
POWER RECLINERS  
By Pride<sup>®</sup> Mobility Products Corp.

THE MOST  
COMFORTABLE  
RECLINER  
IN THE WORLD!

## THE 4-MOTOR RECLINER


WITH  
COUNTLESS  
POSITION  
COMBINATIONS






Shown in Merino  
Steel Fabric

# PERFECTA

Retreat from the bustle of the day and embrace quiet comfort and stillness. Make your home a true oasis with the Perfecta Collection from **VivaLift!**<sup>®</sup> Power Recliners. With one press of a button, the full-width power headrest delivers maximum support to your head, neck and shoulders. Smooth Merino fabrics and a button-back design offer a contemporary look.









**Pride Mobility**<sup>®</sup> (US) 800-800-8586 | (CAN) 888-570-1113  
pridemobility.com | prideliftchairs.com |    





# PERFECTA

COLLECTION	PERFECTA v.1
POSITION	Infinite Lay Flat Relaxed Reclined Upward Lift
WEIGHT CAPACITY	375 lbs.
SEAT-TO-FLOOR	19"
SEAT DEPTH	20"
SEAT WIDTH	22"
TOP OF BACK TO SEAT	25"
FULLY PADDED CHAISE	Yes
FOOTREST EXTENSION	Standard
SEATING TYPE	Total Comfort
CHAIR HEIGHT	43"
DISTANCE FROM WALL	15"
OVERALL WIDTH	31"
CHAIR WEIGHT	142 lbs.
WARRANTY	7-Yr. Pro-rated

## FUNCTIONS AND FEATURES:

-  POWER HEADREST
-  POWER LUMBAR
-  USB REMOTE
-  LITHIUM BATTERY BACKUP
-  FOOTREST EXTENSION
-  POWER BACKREST

## POSITION GUIDE:

-  RELAXED
-  RECLINED
-  UPWARD LIFT
-  INFINITE LAY FLAT

## STANDARD FABRIC: MERINO



FAWN

STEEL

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

 = Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.



**Pride Mobility**® (US) 800-800-8586 | (CAN) 888-570-1113

pridemobility.com | prideliftchairs.com |    