

# Protekt® Aire 1500

Adjustable Alternating Pressure Overlay



HCPCS  
E0181

## FEATURES AND BENEFITS:

- Manual air flow adjustment knob.
- Alternating Cycle Time: 10 minutes.
- Safety end flaps secures the pad to the mattress to prevent slipping.
- Lighted On/Off switch.
- Hanging bed hooks.
- Medical grade fire retardant PVC.
- Thickness: 11.5 gauge.
- Size of pad: 35.5"W x 78.5"L.
- Weight Capacity: 300 lbs.
- Warranty: 18 Month Non-Prorated.

80010	Protekt® Aire 1500 System
80012	Bubble Pad Replacement with Flaps

# Protekt® Aire 1000

Alternating Pressure Overlay



HCPCS  
E0181

## FEATURES AND BENEFITS:

- Economical pressure redistribution surface.
- Alternating Cycle Time: 10 minutes.
- Lighted On/Off switch.
- Hanging bed hooks.
- Medical grade fire retardant PVC.
- Size of pad: 35.5"W x 78.5"L.
- Weight Capacity: 300 lbs.
- Warranty: 18 Month Non-Prorated.

80005	Protekt® Aire 1000 System
80011	Bubble Pad Replacement

# Protekt® Aire 2000

Low Air Loss/Alternating Pressure Overlay System



HCPCS  
E0372

## FEATURES AND BENEFITS:

### PROTEKT® AIRE 2000 PUMP

- Low air loss and alternating pressure therapies in one system.
- Air flow dial to obtain custom patient comfort and pressure levels.
- Preset 10 minute alternating cycle.
- Static function for safe and easy patient transfer.
- Low pressure alarm.
- Adjustable hanging brackets.
- UL approved power cord.

### MATTRESS & COVER

- 18 laser holes keeps the patient cool and dry.
- Removable, water resistant, non-shear, anti-microbial nylon cover.
- Transport connector locks in air for easy patient transfers.
- CPR valve provides rapid deflation of mattress.
- Kink resistant tubing.
- Optional raised side rails available.
- Weight Capacity: 300 lbs.
- Warranty: 18 Month Non-Prorated.

80020	Protekt® Aire 2000 36"x80"x5"
80020RR	Protekt® Aire 2000 with Raised Rails