

# Protekt® Aire 6500

**Low Air Loss/Alternating  
Pressure Mattress System  
660 LBS. Weight Capacity**



HCPCS  
E0277

## FEATURES AND BENEFITS:

- Dual compressor with 18 LPM (liter per minute).
- Alternating pressure cycles of 10, 15, 20 & 25 allows clinicians to customize pressure redistribution and skin surface interface.
- Cell-on-cell mattress design prevents patient from bottoming out in case of power failure.
- 10" air cells are constructed of Nylon/PU for added comfort and support.
- Low air loss function controls and manages moisture to prevent skin maceration.
- 3 alarms; low pressure, power failure, and alternating failure alerts clinicians when preset pressure levels decrease.
- Seat inflate increases air flow by 20% to middle bladders to prevent bottoming out.
- Static mode allows for easy transfer in and out of bed.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Lock out feature prevents unauthorized changes to the pump settings.
- Nylon quilted top cover is fluid resistant, breathable, vapor permeable and fire retardant.
- CPR pull tab quickly releases air from mattress.
- Weight Capacity: 660 lbs.
- Warranty: 2 Year Non-Prorated.

86500	Protekt® Aire 6500 36"x80"x10"
86500RR	Protekt® Aire 6500 with Raised Rails

# Protekt® Aire 7000

**Lateral Rotation/Low Air Loss/  
Alternating Pressure and Pulsation  
Mattress System**



## FEATURES AND BENEFITS:

- Lateral rotation offers three turning modes with adjustable angles and turning cycle times from 5-95 minutes – Left Turn, Right Turn and Full Turn.
- TrueCare™ Alternating enhances the comfort for the patient by fully re-inflating all sets of air cells before the next set of air cells deflates.
- Innovative iLAL™ (independent Low Air Loss) effectively reduces the moisture and skin heat on the support surface by periodically providing a high volume of air loss relief.
- Selectable 10, 15, 20, 25, 30 minute alternating cycle times with LED display and static mode for easy patient transfer.
- Constant low pressure (CLP) mode puts the mattress in a low pressure static surface for patients who cannot tolerate a dynamic surface but still require a pressure redistribution surface.
- Pulsate mode with adjustable cycle time provides a periodic pressure changing static surface.
- Cell-on-cell mattress design prevents patient from bottoming out in case of power failure.
- Auto-firm mode provides quick set up and ease of ingress and egress from bed.
- Audiovisual alarms: power failure, low pressure and service alarm with error code.
- Nylon quilted top cover is fluid resistant, breathable, vapor permeable and fire retardant.
- Two-stage fowler detector automatically stops turning at 20 degrees and increases the pressure output at 45 degrees to prevent "bottoming out."
- Auto-lockout ensures patient safety by locking the panel to prevent accidental setting changes.
- Weight Capacity: 660 lbs.
- Warranty: 2 Year Non-Prorated.

80070	Protekt® Aire 7000 36"x80"x8"
80070-42	Protekt® Aire 7000 42"x80"x8"