

Protekt® Aire 8600AB



**Low Air Loss/
Alternating
Pressure
Bariatric
Mattress System
with High Flow
Pump, Raised
Air Bolsters and
Cell-On-Cell
Support Base**

FEATURES AND BENEFITS:

- High air flow 16 liter pump offers a 660 lb weight capacity.
- Raised side perimeter air bolsters mitigates accidental patient roll-outs, and reduces the risk of entrapment.
- Center cut-out in the air bolsters enables safe ingress and egress from the mattress.
- Alternating pressure cycles of 10, 15, 20 & 25 minutes allows clinicians to customize pressure redistribution and skin surface interface.
- Cell-on-cell mattress design prevents patient from bottoming out in case of power failure.
- 8" air cells are constructed of Nylon/PU for added comfort, support and durability.
- 18 laser air holes to keep the patient cool and dry.
- Audible low pressure alarms alerts clinicians when preset pressure levels decrease.
- Static mode allows for easy transfer in and out of bed.
- Seat inflate increases air flow by 20% to middle bladders to prevent bottoming out.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Lock out key completely locks the control panel, to eliminate tampering of the patient settings.
- CPR pull tab quickly releases air from mattress.

SPECIFICATIONS

- Air Flow: 16 liters per minute
- UL approved power cord.
- Mattress Dimensions: 42"x80"x8" & 48"x80"x8"
- Weight Capacity: 660 lbs.
- Warranty: 18 Month Non-Prorated.

86080AB-42	Protekt® Aire 8600AB w/ Raised Air Bolsters 42"
86080AB-48	Protekt® Aire 8600AB w/ Raised Air Bolsters 48"

Protekt® Aire 9900

**True Low Air Loss Mattress System with
Alternating Pressure and Pulsation**



FEATURES AND BENEFITS:

- 1200 liters per minute true low air loss.
- Provides effective true low air loss with alternating pressure and pulsation.
- State of the art micro-controller technology, monitors blower speed to set pressures.
- User-friendly control panel with dome keys, LED display for simple user interface. Displays patient comfort pressure level by weight.
- Alternating pressure cycles of 5, 10, 15, 20, 25 & 30 minutes allows clinicians to customize pressure redistribution and skin surface interface.
- Cell-on-cell mattress design prevents patient from bottoming out in case of power failure.
- Audio/visual alarms in the event of power outage or low pressures. Internal memory always retains current settings during power outage, will resume desired operation upon power return.
- Lock out key completely locks the control panel, to eliminate tampering of the patient settings, including power switch.
- Optional air bolster helps mitigate accidental roll outs.
- Optional raised side rails available on all sizes.
- Mattress dimensions: 36"x80"x10".
- Available in 42", 48", 54" & 60"
- Weight capacity: 1,000 lbs.
- Warranty: 18 Month Non-Prorated.

81090-36	Protekt® Aire 9900
81090-36RR	Protekt® Aire 9900 with Raised Rails
81090-36AB	Protekt® Aire 9900 with Air Bolsters
81090-42	Protekt® Aire 9900 42"
81090-48	Protekt® Aire 9900 48"
81090-54	Protekt® Aire 9900 54"
81090-60	Protekt® Aire 9900 60"