

# *Protekt™ Aire 9900*

**True Low Air Loss Mattress System  
w/Alternating Pressure & Pulsation**



 **Proactive**  
medical products

270 Washington Street Mount Vernon, NY 10553

**Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005**

[www.proactivemedical.com](http://www.proactivemedical.com) • [sales@proactivemedical.com](mailto:sales@proactivemedical.com)

The Protekt Aire 9900 is the most effective means available of multi-zone pressure therapy, continuous true low air loss relief and pulsation all in one dynamic air system

## FEATURES

- State of the art micro-controller technology, monitors blower speed to set pressures. Software flexibility allows customer changes (custom pressures, times, and functions).
- User-friendly control panel with dome keys, LED display for simple user interface. Displays patient comfort pressure level from Soft 1 to Firm 10 ( $8 \pm 4$  to  $32 \pm 6$  mmHg).
- High air flow blower with 1000 liters per minute.
- Provides effective true low air loss with or without alternating pressure or pulsation.
- Separate dome key for bariatric patients.
- Lock out key completely locks the control panel, to eliminate tampering of the patient settings, including power switch.
- Audiovisual alarms in the event of power outage or low pressures (hose disconnect). Alarm silence key to mute alarm. Internal memory always retains current settings during power outage, will resume desired operation upon power return.
- Therapeutic support surface consists of superior quality, urethane coated nylon air cushions to provide effective alternating pressure therapy and pulsation with true low air loss.
- **OPTIONAL SIDE BOLSTERS** constant air pressure side bolsters help contain the patient on the mattress. Air bolsters can be deflated for easy patient access into and out of the mattress, and also for easy patient procedures.
- Low friction and low shear force producing, urethane coated breathable nylon top sheet, which is liquid proof and highly vapor permeable.
- The air cells are connected to the manifold using quick connectors and fastened to the base. The top sheet is strapped to the base.
- Single rugged, quick disconnect gang connector for rapid attach and CPR deflate functions and transport.
- Complete support surface is made of superior quality materials and is modular in construction. All of the components such as manifold, tubing assembly, air cushions, top sheet, and foam base are interchangeable, and can be easily cleaned or detached for laundry.
- **Weight capacity: 1,000 lbs.**
- **18 Month Warranty.**

81090-P	Protekt Aire 9900 Pump only	Low Air Loss/Alternating Pressure/Pulsation Pump
81090-36	Protekt Aire 9900 Standard 36"	Low Air Loss/Alternating Pressure/Pulsation System 36"
81090-42	Protekt Aire 9900 42" Wide	Low Air Loss/Alternating Pressure/Pulsation System 42"
81090-48	Protekt Aire 9900 48" Wide	Low Air Loss/Alternating Pressure/Pulsation System 48"
81090-54	Protekt Aire 9900 54" Wide	Wide Low Air Loss Alternating Pressure/Pulsation System 54"
81090-60	Protekt Aire 9900 60" Wide	Wide Low Air Loss Alternating Pressure/Pulsation System 60"



270 Washington Street Mount Vernon, NY 10553

**Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005**

**www.proactivemedical.com • sales@proactivemedical.com**