

# PROTEKT® LIFT

## 30600-PLE | FULL BODY PATIENT LIFT

Designed to provide the perfect combination of patient safety, weight capacity and maneuverability. This lift is manufactured to meet the special needs of the bariatric patient population in hospitals and nursing homes. Caregivers appreciate the ease of operation and patients are reassured by the stability and overall comfort provided.

### SPECIFICATION FOR 30600-PLE

weight capacity	600lb
spreader bar highest position	72"
spreader bar max. reach	64"
spreader bar lower position	28"
max. reach at 600mm	23"
max. reach from base	26"
reach from base with legs spread to 700mm	21"
max. height of CSP.	81"
min. height of CSP.	37"
hoist range	44"
max. external width	38"
min. internal width	21"

**WARRANTY:**  
Electronics 2 years    Frame 3 years

**SCALE OPTION:**  
Digital scale available on all PLE lifters. Comes with all necessary hardware.



- Bariatric lift with extra large lift frame provides outstanding integration with large patients.
- 600 pound weight capacity with manual emergency lowering.
- 6 sling hooks allow for multiple sling options and a wider range of transfers.
- Battery lasts twice as long as competing brands for less downtime (24 volt).
- Ergonomically designed foot pedal offers hands free opening of base.
- Clip free sling hookups improve safety and reduce downtime.
- Full parts breakdowns and modular design on lift and electronics.
- Advanced electronics with auxiliary up/down lift controls.