PROTEKT® 500 LIFT

33500 | POWER PATIENT LIFT

The high-performance Protekt® 500 Lift is the perfect blend of form and function for safe and easy patient transfers. This electric-powered, high quality patient transfer system provides an efficient method to safely move patients, preventing patient and caregiver injuries. Designed with the facility in mind, the lift is capable of floor lifts and fits underneath most beds in the lowest position, making it the perfect solution for transferring patients in small rooms with low beds.

- 6 point spreader bar allows for multiple sling options, safer patient handling and a wider range of transfers.
- Capable of lifting from the floor.
- Ergonomically designed foot pedal offers hands free opening of base.
- Manual emergency lowering provides safe response to power loss.
- 4" rear casters with brake & 3" dual front casters.
- Clip free sling hookups improve safety and reduce downtime.
- 500 lb. weight capacity.
- Modular design allows easy parts replacement.
- Advanced electronics with auxiliary up/down lift controls.



(Optional Scale Item # 30300)

WARRANTY: Electronics 2 years Frame 3 years Battery 1 year

SPECIFICATION FOR ITEM # 33500

Weight Capacity	500 lbs.
Unit Weight	109 lbs.
Base Length	47.3"
Base Ext. Width - Closed	26.4"
Base Int. Width - Opened	39.3"
Mast Height	51.3"
Front Caster Height	4.7"
Lifting Height - Lowest	21.9"
Lifting Height - Highest	71.3"
Base Opened Turning Diameter	62.8"
Base Closed Turning Diameter	58.3"

Battery	LP-TC12-BPAK
Cl (DL)	1 D TD7C 4 DD
Charger (Plug)	LP-TP7C-ADP
Charging Cradle (No plug)	LP-TBC2-CRDL
Hand Pendant	LP-TH10-PEN





270 Washington Street, Mount Vernon, NY 10553

Tel: 855-237-7622 • Fax: 845-205-6005

www.proactivemedical.com