

Chariot III Wheelchair



The **Chariot III Lightweight Wheelchair** offers all the benefits of a standard K3 wheelchair with the added feature of removable flip back height adjustable arms and 2" seat extension. Weighing less than 36 pounds (excluding front riggings) the Chariot III also comes standard with tool-free adjustable front riggings and easy push-lock wheel locks. Available with swing-away footrests or elevating legrests.



Removable flip-back armrests are height adjustable



2" Seat extension



Swing-away footrests with heel loops

FEATURES

- Removable flip back height adjustable arms allow for easy egress or ingress
- Dual axle provides easy transition of seat height to hemi-level
- Includes tool-free swing-away footrests with heel loops or elevating legrests with composite footplates
- Elevating legrests are equipped with a pneumatic piston for infinite elevating positions
- Easy-to-clean, comfortable "Nylon" upholstery resists mildew and bacteria
- Cross-brace, powder-coated silver vein steel frame
- 2-inch removable tool-free seat extension
- Easy push-to-lock wheel locks
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability
- Composite handrims with 24" solid urethane tires
- 8" x 1" solid rubber front casters adjustable in 3 positions
- Chart pocket on back for additional convenience
- HCPCS Code: K0003 / K0195 OR E0990 (ELR only), E0973 (adjustable height and detachable arms)

SPECIFICATIONS

Seat Width:	16"/18"/20"
Seat Depth:	16" (18" w/Seat Extension)
Seat-to-Floor Height:	19.5"/Hemi 17.5"
Back Height:	18"
Back Style:	Fixed
Overall Width w/Arms:	24"/26"/28"
Weight Capacity:	300/300/350 lbs.
Product Weight:	36/37/38 lbs.
Shipping Weight:	42/43/44 lbs.
Warranty:	Limited 3 Year

Item #	Description
WCK316AHDAELR	Chariot III 16" K3 Wheelchair w/Removable Adj. Height Flip-Back Desk Arms, ELR, Seat Ext - 300 lb. Wt-Cap.
WCK316AHDASF	Chariot III 16" K3 Wheelchair w/Removable Adj. Height Flip-Back Desk Arms, SWFR, Seat Ext - 300 lb. Wt-Cap.
WCK316AHFAELR	Chariot III 16" K3 Wheelchair w/Removable Adj. Height Flip-Back Full Arms, ELR, Seat Ext - 300 lb. Wt-Cap.
WCK316AHFASF	Chariot III 16" K3 Wheelchair w/Removable Adj. Height Flip-Back Full Arms, SWFR, Seat Ext - 300 lb. Wt-Cap.
WCK318AHDAELR	Chariot III 18" K3 Wheelchair w/Removable Adj. Height Flip-Back Desk Arms, ELR, Seat Ext - 300 lb. Wt-Cap.
WCK318AHDASF	Chariot III 18" K3 Wheelchair w/Removable Adj. Height Flip-Back Desk Arms, SWFR, Seat Ext - 300 lb. Wt-Cap.
WCK318AHFAELR	Chariot III 18" K3 Wheelchair w/Removable Adj. Height Flip-Back Full Arms, ELR, Seat Ext - 300 lb. Wt-Cap.
WCK318AHFASF	Chariot III 18" K3 Wheelchair w/Removable Adj. Height Flip-Back Full Arms, SWFR, Seat Ext - 300 lb. Wt-Cap.
WCK320AHDAELR	Chariot III 20" K3 Wheelchair w/Removable Adj. Height Flip-Back Desk Arms, ELR, Seat Ext - 350 lb. Wt-Cap.
WCK320AHDASF	Chariot III 20" K3 Wheelchair w/Removable Adj. Height Flip-Back Desk Arms, SWFR, Seat Ext - 350 lb. Wt-Cap.
WCK320AHFAELR	Chariot III 20" K3 Wheelchair w/Removable Adj. Height Flip-Back Full Arms, ELR, Seat Ext - 350 lb. Wt-Cap.
WCK320AHFASF	Chariot III 20" K3 Wheelchair w/Removable Adj. Height Flip-Back Full Arms, SWFR, Seat Ext - 350 lb. Wt-Cap.

270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) * Fax: 845-205-6005

www.proactivemedical.com * sales@proactivemedical.com