

Protekt™ Aire 8000BA

Alternating Pressure Bariatric Mattress System



 **Proactive**
medical products

Medicare
Code
EO277

270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com

The Protekt™ Aire 8000BA Bariatric Mattress Series is the first choice by clinicians to treat and prevent bariatric patients from pressure ulcers.

FEATURES

- Alternating pressure cycles of 10, 15, 20 & 25 allows clinicians to customize pressure redistribution and skin surface interface.
- Low air loss function controls and manages moisture to prevent skin maceration.
- 3 alarms; low pressure, power failure, and alternating failure alerts clinicians when preset pressure levels decrease.
- Seat inflate increases air flow by 20% to middle bladders to prevent bottoming out.
- Static mode allows for easy transfer in and out of bed.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Lock out feature prevents unauthorized changes to the pump settings.
- 10" air cells are constructed of Nylon/PU for added comfort and support.
- Quilted top cover for additional patient comfort.
- Audible low pressure alarms alerts clinicians when preset pressure levels decrease.
- CPR pull tab quickly releases air from mattress.

SPECIFICATIONS

- Air Volume: 15 liters per minute
- Cycle Time: 10, 15, 20 and 25 minutes, selectable.
- Pump Weight: 8.8 lbs.
- Pump Dimensions: 13.4"x4.9"x9.7".
- Mattress Dimensions: 42"x80"x10", 48"x80"x10".
- Cell Depth: 10".
- UL/CSA approved power cord.
- Power: 120 VAC.
- Input Frequency: 60Hz.
- Pressure Range: 20-70 mm Hg.
- Weight Capacity: 650 lbs.

80080	Protekt™ Aire 8000BA-42	Alternating Pressure Mattress System-42"
80081	Protekt™ Aire 8000BA-42 Pump	Pump only for Protekt Aire 8000BA-42 & 8000BA48
80082	Protekt™ Aire 8000BA-42 Mattress	Mattress only for Protekt™ Aire 8000BA-42
80083	Protekt™ Aire 8000BA-42 Cover	Cover only for Protekt™ Aire 8000BA-42
80084	Protekt™ Aire 8000BA-42 Bladder	Individual Bladder for Protekt™ Aire 8000BA-42
80085	Protekt™ Aire 8000BA-48	Alternating Pressure Mattress System-48"
80087	Protekt™ Aire 8000BA-48 Mattress	Mattress only for Protekt™ Aire 8000BA-48
80088	Protekt™ Aire 8000BA-48 Cover	Cover only for Protekt™ Aire 8000BA-48
80089	Protekt™ Aire 8000BA-48 Bladder	Individual Bladder for Protekt™ Aire 8000BA-48



270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com