

PROTEKT® TAKE-A-LONG

33400P | FOLDING ELECTRIC PATIENT LIFT

The Protekt® Take-A-Long Folding Electric Patient Lift represents a new standard in full body lifts for individuals who travel or have storage issues in their home. The small foot print is ideal for home use. It easily folds up to fit in most automobile trunks for transport or travel.

- **400 lb. weight capacity.**
- **2 point spreader bar.**
- **Portable and foldable for easy storage and transportation.**
- **Lightweight and easy to move (84 lbs).**
- **Dual rear locking casters.**
- **Tool-Free assembly and disassembly.**
- **Battery powered with a high-capacity battery pack.**
- **Low battery and overload warning alarm.**
- **Built in battery gauge.**
- **Emergency stop button.**
- **Optional spare battery and wall charger available.**
- **Optional digital scale.**

SPECIFICATION FOR ITEM # 33400P

Maximum Capacity	400 lbs.
Base Legs Length	39.4"
Base Closed External Width	25.6"
Base Open External Width	38.2"
Mast Height	50.3"
Base Height	4.3"
Min. Lifting Height	10.6"
Max. Lifting Height	65"
Net Weight	84 lbs.
Shipping Weight	99 lbs.

WARRANTY:

Electronics 2 years
Frame 3 years

