



Bath Stool

The **ProBasics Bath Stool** offers a space saving alternative to standard bath benches while providing the safety patients are looking for. The aluminum frame is lightweight, corrosion-proof and easily disassembled for quick storage. The legs are height adjustable in one-inch increments and feature slip-resistant tips for stability.

FEATURES

- Ideal for narrow tubs and showers
- Texture composite seat is impact resistant
- Manufactured with one-inch aluminum tubing to ensure lightweight stability
- Tool-free, easy to install legs
- Seat height adjusts in one-inch increments
- Color: white
- HCPCS Code: E0242*

SPECIFICATIONS

Seat Diameter:	12.5"
Overall Footprint:	17.5" L x 17.5" W
Adjustable Seat Height:	13.5" – 20.5"
Weight Capacity:	250 lb
Carton Dimensions:	16" x 13.5" x 13.5"
Shipping Weight:	14.5 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSST	Height Adjustable Bath Stool	815067072040	4/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.