

DME



## **Economy Two-Button Folding Steel Walker**

ProBasics Economy Two-Button Folding Steel Walker provides maximum strength in a durable design, all backed by a limited lifetime warranty. Ample height adjustment options coupled with a two-button folding mechanism make it easy to use and accommodating to the most demanding patient needs. Available with or without 5" wheels and in junior or adult sizes. Supports patient weights up to 350 pounds.





Height Adjustable

Push Button Provides Audible Click When Locked

## **SPECIFICATIONS**

SPECIFICATIONS			
Width Inside Hand Grips:	17 "		
Width Inside Base Legs:	20"		
Depth Opened at Base	17"		
Depth Folded:	4.5"		
Wheel Size:	5"		
Patient Height Range:	Junior: 4'4" – 5'7" Adult: 5'2" – 6'4"		
Weight Capacity:	350 lb		
Product Weight:	WKSJN2B: 6.5 lb WKSAN2B: 7 lb WKSJW2B: 7.2 lb WKSAW2B: 7.6 lb		
Carton Dimensions:	WKSJN2B: 21.7" x 9.6" x 26" WKSAN2B: 21.7" x 9.6" x 33" WKSJW2B: 21.3" x 10.6" x 24" WKSAW2B: 21.9" x 10.6" x 31"		
Shipping Weight:	WKSJN2B: 28.8 lb WKSAN2B: 32.6 lb WKSJW2B: 32.5 lb WKSAW2B: 36.2 lb		
Warranty:	Limited lifetime on frame One year on all parts		

ltem #	Description	UPC	UOM
WKSJN2B	ProBasics Economy Two-Button Steel Walker without Wheels, Junior	815067071067	4/cs
WKSAN2B	ProBasics Economy Two-Button Steel Walker without Wheels, Adult	815067071050	4/cs
WKSJW2B	ProBasics Economy Two-Button Steel Walker with 5" Wheels, Junior	815067071043	4/cs
WKSAW2B	ProBasics Economy Two-Button Steel Walker with 5" Wheels, Adult	815067070985	4/cs

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



## **FEATURES**

- 1-inch steel frame construction provides superior strength
- Easy-to-use push button folding mechanism provides audible click when locked
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Front cross brace near top of walker allows for full stride
- Side bar adjusts with height to provide additional support
- Adjustable height settings
- Vinyl contoured hand grip for added comfort and long lasting wear
- Rear glide cap allows for walker to slide easily and smoothly over most surfaces (come standard on models with 5" wheels)
- Folds flat for transport and storage
- HCPCS\* Code: E0135 (without wheels) E0143 (with wheels)