



Three-in-One Folding Commode

The **ProBasics Three-in-One Folding Commode** has three versatile functions; it can be used as a toilet safety frame, a commode, and a raised toilet seat. This full-size steel commode folds compactly for storage and makes delivery easy. Clip-on seat can be removed for easy cleaning. Comes complete with a commode bucket, splash guard and cover.

FEATURES

- Durable steel-welded frame is attractive, strong and easy-to-clean
- Comfortable elongated seat for ease of use
- Easily folds for transportation and storage
- Clip-on seat can be removed for easy cleaning
- Use as a bedside commode, toilet safety frame or raised toilet seat
- Seat height readily adjusts in one-inch increments
- Comes complete with seat, lid, commode bucket with handle, cover and splash guard
- Tool free assembly
- Color: gray frame, gray accessories
- HCPCS Code: E0163*

SPECIFICATIONS

Width between arms:	17.75"
Outside width:	22"
Seat dimensions:	13.5"W x 16.5"D
Seat opening:	8"W x 10"D
Adjustable seat height:	16" – 22"
Weight Capacity:	350 lb
Carton Dimensions:	30" x 23" x 15"
Shipping Weight:	62.5 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSFC	ProBasics Three-in-One Folding Commode with Elongated Seat	815067071845	4/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.