



Low Profile Rollator

The **ProBasics Low Profile Rollator** features a seat that is 3" lower than a standard rollator making it the perfect model for individuals of shorter stature. Constructed of lightweight aluminum, it offers great maneuverability and stability. Equipped with large 6" casters, easy-to-operate ergonomic hand brakes, a water-resistant nylon pouch and a fold-strap for easy storage and transport. Designed to meet all your patients' fundamental needs and height adjustable for a custom fit. Supports patient weights up to 300 pounds.



Water-Resistant
Nylon Pouch



Easy-to-Operate
Ergonomic Hand Brakes



6-inch Wheels

FEATURES

- Lightweight aluminum black frame
- Seat is 3" lower than a standard rollator to accommodate individuals of shorter stature
- Two-button removable curved padded backrest
- Comfortable, flip-up padded seat
- Easy-to-operate ergonomic hand brakes
- Height-adjustable for a custom fit
- Water-resistant nylon pouch attaches under seat for added privacy and security of personal items
- Equipped with a folding mechanism for easy storage and transport
- Non-marring black casters with soft grip tires, ideal for both indoor and outdoor use
- HCPCS* Code: E0143
E0156 (seat)

SPECIFICATIONS

Overall Width:	23.3"
Width b/w Hand Grips:	18"
Seat Dimensions:	14" (W) x 13" (D)
Seat-to-Floor Height:	19"
Weight Capacity:	300 lb
Wheel Size:	6"
Product Weight:	12.5 lb
Carton Dimensions:	23.8" x 7.9" x 31.3"
Shipping Weight:	15 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	Color	UPC	UOM
RLAH6BK	ProBasics Low Profile Rollator	Black	815067071708	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.